



Product insights you can trust.

BRILLIANT EverGlow

EDITOR'S CHOICE



RATING SYSTEM: Excellent + + + + + Very Good + + + + Good + + +

Coltène/Whaledent
www.everglow.coltene.com



96% overall rating

Description

BRILLIANT EverGlow is a highly esthetic and user-friendly universal submicron hybrid composite. Its filler technology is designed for gloss properties such as a high polishability and long-lasting high luster. It also has convenient sculptability without slumping and good wettability for easy handling. The **BRILLIANT EverGlow** shade system has three translucency levels and exceptional blending properties allowing for highly esthetic and flexible single-shade and multi-shade restorations. The medium-translucent Universal shades provide an exceptional blending effect, making them ideal for seamless, single-shade restorations. Duo Shades, which match two VITA shades in one (e.g. A2/B2), help reduce decision time and the amount of material in storage. For use in multi-shade restorations, the Universal shades can be combined with opaque masses and translucent enamel masses. Each **BRILLIANT EverGlow** shade is available separately as a syringe (3 g) or in boxes with 20 individual tips (0.2 g each). Kits containing 4 or 9 shades are also provided and include a Shade Guide with original material for free.

Special Features

- Excellent sculptability, smooth consistency and non-slumping
- High polishability and gloss retention
- Adheres to the bonded tooth surfaces rather than sticking to the instrument
- Esthetic single-shade restorations
- Versatile shade system with three translucencies

Indications

- All anterior and posterior direct composite single and multi-shade restorations
- Luting and repairing of ceramic and composite restorations

"SUPERB LIFE-LIKE ESTHETICS AND LUSTER."

Consultants' Comments

- "Excellent consistency and handling. It is easy to dispense, flows well and doesn't slump."
- "Packs well, good marginal adaptation and can be sculpted."
- "Great chameleon-like properties and one of the easiest composites to finish that I've tried."
- "Great polishability!"
- "Blends well. The blended shade is brilliant for that one patient who's in between shades."
- "Love all the shade options."
- "A little hard to dispense and adapt well."
- "Layers were challenging without having visible lines/seams."



Clinical Tips


- The opaque shades match/adapt well with porcelain repairs and endodontic access opening restorations.
- For shade matching, extrude a little bit on the tooth and light cure it to check.
- Great for teeth with minimum remaining coronal structure.

Evaluation Highlights


BRILLIANT EverGlow was evaluated by 33 consultants, and used 771 times.

- Versatile system with multiple translucency levels
- Easy to handle and place
- Good adaptability
- Esthetic and high luster
- Unique hard plastic syringe sleeve to avoid cross contamination and damage to label when disinfecting. Syringe sleeve may be wiped down with a surface disinfectant

Key Features:



Compared to Competitive Products:



Percentage of Consultants Who Would:

