## 1. Definition

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 **COLTENE**

### 1. Definition

DIATECH Z-Rex are diamond burs for high performance ceramics as for example zirconia and aluminum oxide.

### 2. Indications

- Endo access
- Crown adjustment
- Crown cutting

### 3. Contraindications

- Metal
- Amalgam

### 4. Safety Information

#### 4.1. WARNING

- Jamming of the instrument can cause breakage. See 5.
- Corrosion or strongly acidic and alkaline disinfectants may adversely affect the instruments. See 7.1 & 7.2
- Dry heat sterilization over 140 °C / 284 °F damages the bur. See 7.4
- Tooth structure, pulp and adjacent fillings may be harmed by heat development during preparation. See 5. & 6.
- If not cleaned and sterilized properly, there is a risk of infection. See 7.
- Only supplied to dentists and dental laboratories or upon their instruction.

#### 4.2. CAUTION

- Avoid simple cold disinfectant solutions. These solutions do not sterilize and may contain corrosive agents that degrade the performance and strength of the instrument.
- Only use a cleaning agent/disinfectant suitable for dental rotary instruments.
- Use a cleaning agent/disinfectant containing a corrosion inhibitor.
- The cleaning agent/disinfectant manufacturer's instructions must be obeyed.
- Use sterilizing devices according to the manufacturers recommended procedure. It is the responsibility of the user to ensure that sterilization is effective.
- It is recommended to use a dental dam during a procedure.
- Only use flawless diamond instruments.
- Only use with sufficient amounts of water cooling (minimum 50ml/min).

### 5. Preparations before use

- △ It is recommended to sterilize diamond instruments prior to the first use (see 7. Sterilization)
- △ Instruments must be visually checked: Instruments which are blunt, damaged, bent or no longer concentric, must be removed immediately and disposed of carefully.
- △ An optimal cutting performance and cutting life will be achieved by choosing the correct speed and pressure (depending on the substrate).

### 6. Recommendations for Use and Safe Operation

The diamond instruments have to be selected (shape, size, type) according to the type of preparation. The ergonomic principles must be

taken into account during the working procedures. The high-speed turbine or contra-angle handpiece as well as rotary instruments must be in perfect working condition. The instruments must be fully inserted and gripped in the high-speed turbine or contra-angle handpiece with utmost care and checked for imbalance and concentricity.

Before touching the preparation site with the instrument, make sure that the bur is rotating at the optimum speed (see packaging) and that the cooling system is working properly. Once the preparation has been completed, remove the rotary instrument from the site and allow it to come to a standstill.

The recommended force applied (reference value 1.5 N) is based upon the diameter of the instrument, the type of handpiece and the substrate to be worked on. Use recommended speeds according to the table under "Recommended speeds". Use lower speed and pressure for finer grit size or for finishing, and a higher speed and pressure for coarse grit sizes and for bulk reduction. Only use with water cooling (minimum 50ml/min).

#### Recommended speeds

Head diameter in 1/10 mm	Speed range (RPM)
012	160 000 – 300 000
014	135 000 – 275 000
016	120 000 – 240 000
018	105 000 – 210 000
023	85 000 – 165 000

### 7. Cleaning after use

**7.1. Manual Cleaning:** Presoak in enzymatic cleaner (e.g. BioSonic UC32) to loosen debris (5 min.). Brush away re-maining debris. Rinse (2 min.), dry with absorbent lint free towel.

**7.2. Ultrasonic Cleaning:** Burs may be ultrasonically cleaned after inserting them in bur holders (or blocks) to prevent damage. A cycle of 10 minutes is recommended, using an enzymatic cleaner (e.g. BioSonic UC40), followed by a cold water rinse (2 min.). Be sure to follow manufacturers' in-structions.

**7.3. Inspection after Cleaning:** Visually inspect bur, with magnification if needed, to ensure that all debris has been removed; if not repeat cleaning process. **WARNING:** Burs cannot be sterilized properly unless they are thoroughly cleaned; free of debris.

#### 7.4. Sterilization


**Autoclaving:** Burs may be sterilized by using a dynamic air removal sterilization cycle. Sterilize in a pouch (e.g. SPSmedical Self-Seal Pouches) at full cycle with a dwell at 132 °C / 270 °F minimum for 3 minutes. Use the sterilization device according to the manufacturer's IFU.

### 8. Storage

Dental burs should be stored in a moisture free environment. Whenever possible, the

instruments should be stored in their original packaging to enable them to be identified and traced.

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COLTENE/WHALE DENT AG