
SciCan HYDRIM® L110W G4 Quick Reference Guide

Operator's Manual should be consulted before using this washer.

BASIC OPERATION

1 Power unit ON.

2 Load Instruments.

3 Close the door.

4 Select cycle and press START.

TOUCHSCREEN

Rinse and Hold Cycle

Prevents soil from drying on instruments. Incomplete loads that will not be washed within the hour. NOT a cleaning cycle.

Regular Cycle

Moderately soiled loose instruments.

Heavy Duty Cycle

Heavily soiled instruments and cassettes.

Door Lock

User/Technical Menu

Network Connection

USB

Water Softener/Detergent Level/Door Status


Enter

Back

Refilling the Water Softener Salt Reservoir

1. Unscrew lid.
2. Fill to the top (max. 1 kg / 2.2 lbs.)
3. Screw lid on tightly.
4. Run a complete cycle with no instruments.

IMPORTANT!

If the salt container is being filled for the first time, it MUST be filled with at least 1 litre of water before adding the salt.

Only tablet salt should be used in the softening unit. Fine grained salt such as table salt should NOT be used as it could clog the unit.

The water softener salt is also available from SciCan. Please contact your local dealer.

Cleaning the Chamber

Pour 1.0L of vinegar or citric acid-based instrument neutralizer into the chamber and follow these steps:


IMPORTANT! Do a Chamber Cleaning cycle every 25 cycles. Failure to clean the chamber can promote corrosion.

Replacing the HIP™ Cleaning Solution

Recommended precautionary measures with regard to safe handling of HIP solution:

- Wear undamaged surgical gloves and eye protection.
- After handling the product, remove and dispose of the gloves and wash hands prior to removing eye protection.

TROUBLESHOOTING

Refer to Operator's Manual for a more comprehensive list. Call for service if suggested recommendations do not resolve the error message.

CF2 Chamber Filling Water is not filling the chamber	<ul style="list-style-type: none"> Check that the water hoses are not kinked. Check that water shut-off valves are opened. Call for Service.
CF15 Water Reservoir Overflow Overflow switch triggered	<ul style="list-style-type: none"> Restart the program. Call for Service.
CF21 Dosing System Detergent is not dispensing	<ul style="list-style-type: none"> Tighten detergent cap. Check for leaks. Call for Service.
CF28 No Water Pressure	<ul style="list-style-type: none"> Check the water shut-off valves are open. Check water hoses are not kinked. Call for Service.

Spare Parts and Accessories

CS-HIPL	HIP Cleaning Solution (2 x 3.8L)
01-112594S	Water Softener Salt
01-108305S	Water Test Kits
WTI-SK	Wash Test Indicator - Starter Kit
WTI-RK	Wash Test Indicator - Refill Kit

Terminal sterilization of instruments **MUST** follow processing in the HYDRIM L110W G4.

FOR DETAILED INFORMATION PLEASE SEE PROVIDED OPERATOR'S MANUAL

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