BioSonic® UC300R

Recessed Ultrasonic Cleaning System

Information for User's

- EN** Information for User's
- FR** Guide de l'utilisateur
- ES** Manual del usuario

An Essential Part of Your Infection Control System.

Congratulations on your purchase of a BioSonic Ultrasonic Cleaner — the finest, most efficient and most popular ultrasonic cleaning system in dentistry.

Developed specifically for the dental industry, BioSonic units feature unique product design and optimize time and safety in the office or laboratory.

BioSonic units are a complete ultrasonic cleaning system, with a full line of cleaning solutions and unit accessories which provide you with the flexibility to efficiently clean instruments, prostheses and other items.

What is Ultrasonic Cleaning?

Ultrasonic Cleaning is a process created by high frequency soundwaves. The soundwaves, enhanced by specially formulated cleaning solutions, create high energy cavitation. During cavitation, millions of tiny bubbles form and then collapse or "implode," releasing enormous amounts of energy and shock waves which scour the surface of instruments, appliances and other apparatus placed in the cleaning solution. This powerful scouring action reaches into minute crevices which manual brushing cannot reach. The combination of energy and specially formulated solutions make ultrasonic cleaning the most effective method for removing gross and microscopic debris.

Intended Use: The Biosonic® UC300R Ultrasonic Cleaner is intended for cleaning medical instruments by emission of high frequency soundwaves.

I. Unpacking

1. Remove your BioSonic unit from the shipping container and inspect the unit for any possible damage which may have occurred during shipping.
2. Check the serial numbers on the electronics enclosure and tank enclosure unit with the number on the shipping carton. If they are not the same, contact your dealer.
3. Fill out the warranty card and mail to Coltène/Whaledent within 10 days of receipt.

II. Installation Instructions

The BioSonic UC300R consists of a tank enclosure, an electronics enclosure, a drain tube, a coaxial cable which is connected between the N female coaxial connectors on the electronics control housing and the tank housing; and a 5 pin male / female connector cable which connects between DIN connectors on the tank enclosure and the electronics enclosure, and a power supply cord.

Tools Required:

- 1 Electric Drill
- 1 Drill bit to drill a hole large enough to accept your reciprocating saw blade
- 1 Reciprocating saw and blade

1 Phillips head screw driver

1 Flat blade screw driver

Tank Enclosure:

Note: This Ultrasonic Cleaner must be installed by an experienced equipment manager or licensed contractor. Be sure to follow local building codes.

Note: Before you begin please read and understand the safety instructions provided by the manufacturer of the tools required for installation.

Note: Cables must be installed in the tank assembly prior to the installing the tank assembly into the countertop. Once installed, care must be taken in lining up and installing the other end of the cables to the component box.

1. Select the desired location for the cleaner. You will need at least 11 inches clear vertical space below the top of the counter to accommodate the tank enclosure and an additional 6 inches at the location of the valve for a total of 17 inches. The countertop must be flat and level. Make sure that you are not selecting a location which will interfere with drawers or vertical partitions between cabinets.
2. Tape the enclosed template down to the counter top in the desired location. Make sure that the "front of the unit" is closest to the user.
3. Drill four (4) 3/8" dia. holes where indicated on the template (holes marked "A").
4. Drill four (4) holes either 3/8", 1/2", 3/4", 1" dia. or larger at the four locations marked "B." Select your hole size based upon the saw you will use in "#5" below.
5. Cut along the "cut line" with a jigsaw or other reciprocating saw.
6. Clean away dirt or sawdust around the opening, and remove scrap wood.
7. Place tank enclosure into opening. Be sure that the printed BioSonic logo on the bezel faces the user.
8. Assemble the bearing plates, lock washers and wing nuts to the tie-down studs (per Detail A).
9. Tighten all the wing nuts a few turns at a time until the bezel just contacts the counter top. Once the bezel has contacted the counter top on four sides, back off each wing nut 1/4 turn.
10. (Optional) Apply a thin bead of clear RTV silicone sealant between the bezel and counter top.
11. Select a location for the electronics enclosure which will be conveniently close to the tank enclosure.

Note: Mount control unit so there is access to fuses. Be sure to leave at least 4 inches between the back and sides of the electronic control housing and cabinet wall. This will allow for adequate ventilation as well as facilitate the insertion and removal of power/RF cables. Note that the side of the unit has ventilation louvers. Position the electronics control housing in such a way that liquids cannot flow into the unit.

12. Four (4) #8 Phillips head wood screws have been included

if you wish to mount the electronics enclosure to one of the interior cabinet walls. Four adhesive-backed rubber feet have been applied to the underside of the electronic control unit (**Fig. 1**).

13. Connect the two cables between the electronics control box and the tank enclosure. The cables are configured such that they cannot be inserted into the wrong receptacle. If you have difficulty plugging in a cable, check that it matches the receptacle.
14. Plug the power cord into the electronics control box, then into a 115V wall outlet.

	Power should not be supplied to the unit until all interconnections between the electronics enclosure, tank enclosure and any accessories are securely connected.
	<p>WARNING! (Failure to comply may result in personnel injury and/or equipment damage.)</p> <p>This unit must be properly grounded. Failure to properly ground unit may result in serious injury.</p>

15. The hose supplied with the UC300R will allow you to dispose of the waste in accordance with your local, state and/or federal regulations.

Auxiliary Digital Remote Timer (optional) Purchased separately (**Fig. 12**).

Cat. No. UC315115

Note: The UC315115 Digital Remote Timer is only to be used in conjunction with the Biosonic UC300R 115V Ultrasonic Cleaner and is not intended to be used with other Ultrasonic Cleaner models. The Digital Remote Timer is designed to be permanently mounted into a cabinet face using the supplied electrical box. Installation procedures are detailed below.

Tools Required:

- Electric Drill
- Phillips-Head Screwdriver
- Reciprocating Saw & Blade
- Small Flat-Blade Screwdriver
- ¼" Socket
- Torque-Wrench Capable of in-lbs

	Unplug the cleaner from the power source before beginning installation.
	<p>WARNING! (Failure to comply may result in personnel injury and/or equipment damage.)</p>

Note: The following work must be performed by a licensed electrician and in accordance with local building codes.

Note: The Digital Remote Timer may only be installed using the supplied Carlon (Thomas & Betts®) model B225R in the manner described in this document.

1. Select a convenient location for the timer no further than 92" from the electronic control housing. Assure the location will not interfere with internal cabinet components or other existing installations.
2. Create a template for the electrical outlet box by tracing its outline onto the cabinet face. Cut out the template using the reciprocating saw.
3. Disconnect the cable from the timer terminal block and relay.
4. Install the strain relief adapter into the outlet box knock-out.
5. Apply the strain relief bushing approximately 3 ½" from the end of the cable jacketing. Guide the timer cable into the outlet box and install the strain relief bushing into the knock-out adapter.
6. Slide the electrical box into the cabinet assuring that it sits flush with the surface. If the cabinet thickness is less than 1-1/4" the box locking tabs may be utilized to secure the box, otherwise, install four (4) #6 x 1" screws through the box face and installation surface.
7. Reconnect the timer wiring as follows (Refer to **figure 11**):
 - Green/Yellow ground conductor to the threaded stud located on top of the terminal block. The ring terminal sits on top of the flat washer and beneath the lock washer. Tighten the nut to 3 in-lbs using a 1/4" socket
 - Connect the brown line conductor to the terminal block opposite the "comb" shaped jumper.
 - Connect the blue neutral conductor to the terminal block opposite the pre-connected blue conductor.
 - Connect the brown/white conductor with the quick-connect terminal to the (N O) "normally open" terminal on the timer relay. The connector should be orientated so that the flat side lines up with the a side of the terminal recess.
8. Place the timer into the electrical box. Screw the provided four (4) #6 x 1-1/4" oval head screws through the timer faceplate and into the corresponding holes in the electrical box. Assure that the faceplate is tight against the surface and does not provide access to the wiring beneath.
9. Plug the timer's 4-pin connector into the corresponding 4-pin receptacle on the electronic control housing (**Fig. 13**). Secure the cable to the electronic control housing by screwing the connector onto the threaded receptacle.
10. Using a small flat-blade screwdriver, slide the red timer switch on the back of the electronics housing to read "ext." This will disable the mechanical timer on the electronics control housing and enable the Digital Remote Timer.
11. Plug the cleaner power cord back into the electrical outlet.

You are now ready to begin using your Digital Remote Timer. For proper operation, refer to the "Instructions for Use" provided with the timer.

III. Operating Instructions

A. General Cleaning

1. Fill the unit with warm water to the fill line indicator (**Fig. 2**).
2. Pour the required amount of ultrasonic cleaning solution using either: BioSonic General Purpose Cleaner, BioSonic Enzymatic Cleaner or BioSonic Germicidal Cleaner into the tank.

 <p>CAUTION!</p>	<p>Never operate the unit without solution. Operating a dry tank may result in permanent damage and will void the warranty. Do not place objects directly on bottom of tank, use accessory baskets, cassette racks or beaker holders.</p>
--	---

3. Plug the line cord into a three-wire grounded AC outlet (or per local regulations).
4. To activate the unit, turn the timer to the desired cleaning time.

The indicator light will glow when the unit is in operation. Unit will automatically shut off at the end of the cycle.

5. Change the solution in the tank daily, or sooner if it appears soiled or discolored. When beakers are used, change the solution in the beaker after each use.
6. Drain unit by opening the valve and allow solution to empty into a waste storage container for disposal in accordance with your local, state and/or federal regulations. Rinse tank and close valve (**Fig. 3**).

B. Degassing

Degassing the ultrasonic solution is an important step. Degassing is the removal of large air bubbles found in all fresh solution, which act as a barrier to efficient cleaning.

Note: Degassing is necessary every time new solution is added.

Degas the solution by activating the unit for 15 minutes before you add any items to be cleaned.

C. Cover

The BioSonic Ultrasonic Cleaner should always be covered when operating to trap airborne contaminants in the unit. The BioSonic UC300R is equipped with a "hinged" cover for easy one hand opening. The cover also conveniently lifts off for cleaning or for use in spaces with limited space above the unit.

D. Instrument Draining

The UC300R is equipped with a bezel surrounding the tank. The bezel drains condensation back into the tank keeping the work area around the unit dry. The bezel also acts as a draining stand for the large basket or cassette racks.

Simply position the corner of the basket or rack on two of the draining tabs at opposite corners of the tank (**Fig. 4**) and allow the solution to drain off instruments.

IV. Cleaning Methods

A. In Tank

Basket Cleaning:

Place parts to be cleaned into the basket. Then place the basket directly into the main tank containing either: BioSonic General Purpose Cleaner, BioSonic Enzymatic Cleaner or BioSonic Germicidal Cleaner. Clean according to instructions on bottle label.

Note: Baskets are furnished with adjustable partitions to keep instrument set-ups separated by user or procedure, eliminating sorting and speeding instrument processing. (**Fig. 5**) For more information, see section "V. Accessory Use Directions".

Cassette Cleaning:

Place filled rack into tank containing BioSonic General Purpose Cleaner, BioSonic Enzymatic Cleaner or BioSonic Germicidal Cleaner.

B. In Beaker

For specialized item cleaning e.g., prostheses:

Place item into a beaker, add desired cleaning solution (e.g., BioSonic Cement Remover) according to the manufacturer's directions. Suspend the beaker in the tank by the use of the beaker holder and place the beaker positioning ring around the beaker so that the bottom 2/3 will be suspended in the main bath.

V. Accessory Use Directions

A. Large Instrument Basket:

BioSonic UC310 Large FingerGuard Instrument Basket cleans large volumes of instruments and increases instrument processing efficiency.

The Large FingerGuard Instrument Basket comes with five (5) adjustable basket section inserts to help keep instruments sorted by user or procedure, eliminating time consuming instrument sorting (**Fig. 5**). Additional basket section inserts are available in packets of two (Cat. No. UC335).

Use Directions:

The Large FingerGuard Instrument Basket may be used with or without adjustable section inserts. Place instruments across ridges in bottom of basket. Place basket in tank and clean for desired time.

After cleaning, rinse basket and instruments thoroughly under running water. After rinsing the large FingerGuard basket may be placed on the draining tabs on the bezel for draining.

Note: Beaker holder may be used in conjunction with the large FingerGuard basket (**Fig. 6**). Remove basket section inserts seated directly under beakers to provide space for beakers to seat fully. Use caution when removing basket if beakers are left in place.

B. Twin FingerGuard Basket:

The UC311 Twin FingerGuard Baskets provide you with the versatility of cleaning either two (2) separate basket loads of instruments or a basket and beaker combination for prostheses cleaning (or other small items). Each Twin FingerGuard Basket comes with two (2) basket section inserts.

Use Directions:

Place Twin FingerGuard Baskets side by side (**Fig. 6 & 7**), or use one basket with beaker holder. Twin FingerGuard Baskets may be used with or without basket section inserts. Place instruments across ridges in bottom of basket. Place basket in tank and clean for desired time.

C. Basket Section Inserts:

The UC335 basket section inserts keep instruments separated by user and procedure, eliminating sorting. This increases efficiency and reduces the potential for accidental sticks and cuts since instruments are pre-sorted and handled less.

Use Directions:

Place locking tabs of adjustable basket section inserts into basket holes (**Fig. 8**). Inserts can be placed at any desired width along the basket. Place as many basket section inserts as required for your cleaning needs. To remove or adjust position, simply squeeze center tabs to unlock section insert from basket holes.

D. Beaker Holder:

Place the UC330 Beaker Holder on bezel. Place beaker positioning rings around the beaker so that at least 2/3 of the beaker is suspended in the main bath. Clean for desired time.

E. Cassette Basket:

The UC325 instrument cassette basket may be used to hold up to 3 Large or 6 small standard Hu-Friedy IMS cassettes.

Use Directions:

Place the filled basket (**Fig 9**) into the tank containing either Biosonic General Purpose Cleaner, Biosonic Enzymatic Cleaner or Biosonic Germicidal Cleaner (**Fig. 10**). Once the cassettes are cleaned, remove the basket from the tank and tilt back over the tank to drain. After cleaning rinse the basket and instruments thoroughly under running water. After rinsing the basket may be placed on the draining tabs on the bezel for proper draining. (See instrument draining section for details.)

Section Insert Support Installation Instructions: Place the section insert supports in your desired location along the raised ridge in the bottom of the UC310 Large FingerGuard Instrument Basket or UC311 Twin FingerGuard Basket. The section insert supports are inserted into the holes along the top of the raised ridge with the tapered end facing upward toward the top of the supports. Position the supports at equal distance from either end of the basket, along the raised ridge, to allow the basket section insert to rest squarely in the basket. Slide basket insert into the slots of the section insert supports.

VI. Environmental Conditions

- Intended for indoor use.
- Maximum altitude 2000 m.

•  40°C (104°F)
5°C (41°F)

•  80% @ 31°C (88°F)
50% @ 40°C (104°F)

- Mains supply voltage fluctuations not to exceed $\pm 10\%$ of the nominal voltage.
- Transient over voltages category II.
- Pollution degree 2.
- Class I equipment.

EMC Class A: Suitable for use in all establishments other than domestic and those directly connected to a low voltage power supply network which supplies buildings used for domestic purposes.

VII. Cautions and Warning

NOTE: Failure to comply with these Caution and Warning statements may result in injury to personnel or equipment damage.

 WARNING! (Failure to comply may result in personnel injury and/or equipment damage.)	This unit must be properly grounded. Failure to properly ground unit may result in serious injury.
	Do not place hand or fingers into ultrasonic tank while unit is in operation.
	Do not use flammable liquids with BioSonic units.
	Keep unit covered to reduce airborne contaminants.
	Unspecified or improper use of this product may impair safety protection.

 CAUTION!	Do not operate the unit without sufficient solution in the tank.
	Do not use water alone as the carrier bath. Water is not a satisfactory coupling agent unless a wetting agent is used. For best results use BioSonic General Purpose Cleaner, BioSonic Enzymatic Cleaner or BioSonic Germicidal Cleaner solutions.
	Do not place objects directly on bottom of the tank. They may block cavitation and cause buildup of heat in the transducer resulting in premature failure. Use only BioSonic accessory baskets, cassette racks or beaker to hold items to be cleaned.
	Ultrasonic cleaning of gemstones, pearls and jewelry is not recommended as they may be damaged by ultrasonics.
	Do not restrict air flow. Provide adequate ventilation.
	Do not subject unit to shock or impact.
	Do not immerse the unit in water.

VIII. Order Information

Unit and accessories

CAT No.	Description	Contents
UC300R	115, Unit (300watts)	1
UC310	Large FingerGuard Instrument Basket	1
UC311	Twin FingerGuard Instrument Basket	1
UC330	Two Beaker Holder	1
UC335	Basket Section Inserts	2
UC340	Replacement Cover	1
UC350	Replacement Drain Filters	3
UC355	Replacement Drain Hose	1
UCB2	Small Parts Basket	1
UC53	600 ml Beaker with Cover	1
UC315	Remote Digital Timer	1
UC325	Cassette Basket	1

IX. Tank Solutions:

CAT No.	Description	Contents
UC30	General Purpose Cleaning Solution Concentrate in MeterDose® Dispensing Bottle (16 oz.)	473 ml
UC32	Enzymatic Ultrasonic Cleaner Concentrate in MeterDose Dispensing Bottle (8 oz.)	236 ml
UC34	Plaster & Stone Remover Powder Concentrate in Jar (1.85 lbs.)	840 g
UC38	Germicidal Ultrasonic Cleaning Concentrate in MeterDose Dispensing Bottle (8 oz.)	236 ml
UC39	Cement Remover Concentrate in MeterDose Dispensing Bottle (16 oz.)	473 ml
UC42	Germicidal Ultrasonic Cleaning Concentrate in MeterDose Dispensing Bottle (16 oz.)	473 ml

X. Unit Specifications

Tank capacity: 3.8 gallons (14.4 liters)

Use capacity: 3.0 gallons (11.4 liters)

Overall dimensions:

Above countertop surface:

17.125 x 12.125 x 2.25 inches (43.5 x 30.8 x 5.4 cm)

Below countertop surface:

14.125 x 10.5 x 14 inches (36.2 x 26.7 x 35.6 cm)

Tank dimensions:

14 x 9 x 6.75 inches (35.6 x 22.9 x 17.2 cm)

Consumption: 300 watts max

Operating frequency: 40 kHz

115V Fuse 7A/250V Slow Blow 3AG

XII. Symbols

Caution, Hot Surface

Warning/Caution: Refer to accompanying documents.

Caution: Risk of Electric Shock

Temperature Limits

Humidity Limits

Proper Disposal Required

IEC: 417-5019 Protective conductor terminal

XIII. Maintenance

A. Cleaning

All components should be cleaned each day by wiping with a soft cloth moistened with a commercially available mineral deposit remover followed by cleaning with isopropyl alcohol or a mild disinfectant.

 CAUTION!	Do not use petroleum based solvents, iodophors or phenolic based products. (Iodophors and phenolics can stain the surface of the unit.) Clean up all liquid spills immediately. Rust deposits on tank can be cleaned/ removed with a commercially available stainless steel cleaner.
---------------------	--

B. Stainless Steel Cleaning

Allowable:

- Wash with a clean cloth or rag soaked with warm water and mild soap or liquid dish detergent.
- Always clean in the direction of the original polish lines.
- Always rinse well with clear water (2 or 3 times) after cleaning. Wipe dry completely.
- Specialized household stainless steel cleaner may also be used.

Do Not Use:

- Steel or stainless steel wool or scrapers to remove stubborn debris.
- Harsh or abrasive cleaners.
- Do not let dirt accumulate.

Avoid:

- Cleaners that contain bleach.
- Products that contain chloride, fluoride, iodide, or bromide.

Chemicals harmful to use in the tank:

Acetophenone	Aluminum Chloride
Aluminum Floride	Aluminum Sulphate
Aluminum Biflouride	Ammonium Chloride
Ammonium Hydroxide	Amyl Chloride
Antimony trichloride	Aqua Regia
Bromine	Calcium Bisulfate
Calcium Bisulfite	Calcium Hypochloride
Chloracetic Acid	Chloric Acid
Chlorine Anhydrous	Copper Chloride
Copper Fluoborate	Ethyl Chloride
Ferric Chloride	Ferrous Chloride
Ferris Sulfate	Fluoboric Acid
Flourine	Hydrobromic Acid
Hydrochloric Acid	Hydrocyanic Acid
Hydroflouric Acid	Hydrofluosilicic Acid
Iodoform	Mercuric Chloride
Muriatic Acid	Phosphoric (crude)
Sodium Hypochlorite	Potassium Chloride
Stannic Chloride	Stannous Chloride
Sulfuric Acid	Zinc Chloride

C. Inspecting the Power Cord

Periodically check the power cord for damage and see that the connectors at both ends of the cord are fully seated.

B. Service

This unit must be serviced by an authorized Coltène/Whaledent Service Center and should always be shipped complete with tank assembly and electrical control unit.

XIV. Warranty Information

Our products are carefully manufactured to meet stringent quality assurance requirements. Our products are manufactured from new parts or new and serviceable used parts. Regardless, our warranty terms apply. This product has been developed specifically for use in dentistry and is intended to be operated only by qualified dental professionals in accordance with the instructions contained in this guide. However, notwithstanding anything contained herein to the contrary, the user shall at all times be solely responsible for determining the suitability of the product for the intended purpose and the method of its use. Any guidance on application technology offered by or on behalf of the manufacturer, whether written, verbal or by demonstration, shall not relieve the dental professional of his/her obligation to control the product and to make all professional judgments regarding its use.

Our products are warranted in accordance with the terms of a written Certificate of Limited Warranty accompanying each product. Except for the warranties specifically set forth in the Certificate of Limited Warranty, Coltène/Whaledent Inc. provides no warranties or guarantees of any kind covering the product, expressed or implied, including, without limitation,

any warranties as to merchantability or fitness for a particular purpose. The purchaser/ user is referred to the Certificate of Limited Warranty for all of the terms, conditions and limitations of the warranty covering this product. This Section of the user manual is not intended to in any way modify or add to the warranty provided in the Certificate of Limited Warranty.

Any claim for damage or breakage to the product in transit should be made to the carrier promptly upon discovery. CW does not warrant the product against shipping damage.

XV. Disposal

The Biosonic® UC300R is an electric appliance requiring proper disposal according to local laws and regulations.

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© 2016
Coltene/Whaledent Inc.
Coltene/Whaledent GmbH + Co. KG