

## GENERAL USAGE

### 1 Opening / Closing the Cassette

**5000 Model**

**2000 Model**

With thumbs on the thumb pads, push down, pulling up with forefingers until lid opens.

### 2 Filling Cassette

**2000 and 5000 Model**

**2000 Model**

**5000 Model**

**Do not overload cassette.**

### 3 Inserting Cassette

**5000 Model**

Hold the cassette handle in one hand and the carry handle in the other. Place the end of the cassette into the unit and drop the carry handle into its closed position.

**2000 and 5000 Model**

Gently push cassette inward until you hear a click. Never force the cassette into the STATIM.

### 4 Cycle Starting Guide

**Unwrapped Cycle (Solid)**  
**Unwrapped Cycle (Hollow)**

Sterilise light loads of solid metal instruments such as hand instruments, speculums and forceps; or hollow instruments such as dental/phaco handpieces. (Press button again to view available cycles)

---

**Wrapped Cycle**

Sterilise solid or hollow metal instruments which have been sealed in autoclave bags.

---

**Rubber and Plastics Cycle**

Sterilise solid or hollow unwrapped instruments.

---

**Air Dry Only Cycle**

This is NOT a sterilisation cycle. Refer to Operator's Manual.

## TOUCH SCREEN MENU OPTIONS

**Network Connectivity** – Indicates you are connected to the network.

**USB Storage Device** – Access last 5 cycle data and last 5 cycle faults. Indicator also confirms proper USB connection.

**Water Level / Quality Status\*\*** – Indicator will notify user if water is of sufficient quality or quantity to proceed with sterilisation.

**Cassette Status\*\*** – Indicator will notify user if cassette has been properly inserted.

**User/Technician Menu** – Access unit information and menus to change unit settings.

**START Button** – Button is enabled once a cycle type has been selected and used to begin a cycle.

**STOP Button** – Button is enabled after a cycle has started and used to interrupt a cycle.

Toggle **UP/DOWN** to view menu options.

**ENTER**  
Push button to view/select menu options.

**BACK**  
Push button to return to previous menu.

\*\* NOTE: Sterilisation cycle will not begin if there is an error

## DRYING PERFORMANCE Four steps to obtaining optimal drying results

### 1 Set Up Waste Bottle

Fill this bottle with water to the minimum level and set it on, or close, to the floor.

### 2 Levelling your Unit

1. → →
2. Scroll to **Level** and select.
3. Adjust the leveler feet to move the bubble. Level will change from red to green.

### 3 Coat Cassette

After cleaning and drying the cassette, treat the interior surfaces with STAT-DRI® Plus.

### 4 Use Drying Plate Rack (5000 only)

For optimal drying performance of wrapped instruments, use the drying plate rack supplied and insert one bag per plate with window side up.

Plate positions for full load (10 bags).

## GENERAL MAINTENANCE To ensure trouble-free performance, follow this preventative maintenance schedule.

| Frequency      | Component              | Maintenance Task   |
|----------------|------------------------|--|
| Daily          | Water Reservoir        | Replace water as needed. For ophthalmic use, drain at the end of every workday, leave empty, and refill at the start of the next workday.  |
|                | Waste Bottle           | Empty the waste bottle every time you refill the water reservoir. Fill the bottle with water, up to the MIN line marking. You may also add some chlorine-free disinfectant.  |
| Weekly         | Cassette               | Wash the interior of the cassette with dishwashing soap or a mild detergent that does not contain chlorine. Scrub the inside with a cleaning pad designed for use with Teflon®-coated surfaces. After removing all traces of the detergent, treat interior surfaces of the cassette with the STAT-DRI® Plus drying agent to enhance drying process. Additional STAT-DRI-Plus can be ordered by requesting 2OZPLUS, 8OZPLUS, or 32OZPLUS. |
|                | Water Reservoir Filter | Check the water reservoir filter every week and clean if necessary. Replace only if necessary.   |
| Every 6 Months | Biological Filter      | Replace every 500 cycles or six months (whichever is first). Call for service if wet.  |
|                | Air Filter             | Replace every 500 cycles or six months (whichever is first).   |
|                | Cassette Seal          | Replace every 500 cycles or six months (whichever is first), or whenever necessary.  |

**IMPORTANT!**

Users of this device must be sure to select the most appropriate sterilisation cycle according to the recommendations of the leading infection control authorities and local regulatory guidelines/recommendations.

