

The Monoblock Technique – a revolution in adhesive dentistry?

PRIV. DOZ. DR. MED. DENT. STEFAN J. PAUL

A composite cement with an integrated bonding system that can also be used as a core build up material has been a long-awaited dream in restorative dentistry.

According to standard practice today, 3 to 4 different materials, which are often from different manufacturers, are required for bonding to dentin and enamel, fabricating composite core build ups, and adhesive cementation. Since modern composite materials in dentistry are still based on methacrylate, combining materials – for example, a bond from one manufacturer with the composite of another – is often not a problem. Nevertheless, it would be desirable to have one integrated system available.

Coltène/Whaledent has recently developed a dual-curing composite material that can be used as a cement as well

as a core build up material (ParaPost ParaCore Automix 5 ml). A chemical curing dentin bonding agent, which is compatible with the material, is also available (ParaBond consists of a Non-Rinse Conditioner and Adhesive A & B, which requires mixing before application; and is why it is defined as a two-step bonding system). ParaBond and ParaCore can be used for: 1. adhesive cementation of a root canal post, 2. fabrication of a core build up, and 3. adhesive cementation of a permanent restoration. Coltène/Whaledent describes this time-saving application as the “Monoblock Technique.” The ParaBond/ParaCore System demonstrated excellent sealing against marginal microleakage, which indicates good to very good clinical viability²⁰.

The Monoblock Technique is particularly suitable when light-transmitting, metal-free root canal posts are used with

endodontically-treated teeth that will be fitted with a crown. Root canal posts provide greater retention of the core build up, and distribute masticatory forces along the interface of the residual tooth structure^{3,10,17}. The use of metal-free root canal posts prevent the greyish translucency at the gingival margin caused by the light reflexion from metal root canal posts.

Publications regarding the use of root canal posts recommend that any further weakening of the residual tooth structure caused by the use of a post should be avoided as much as possible. Root canal posts are primarily indicated whenever there is over 50% loss of tooth structure. The smallest diameter of root canal post should always be used to ensure that there is no overloading of the abutment tooth by the final restoration^{4,5,13,19}. From an esthetic point of view, prefer-



Fig. 1: Preoperative clinical situation



Fig. 2: Sufficient root canal filling



Fig. 3: Shade selection by the dental technician



Fig. 4: Preoperative clinical situation with a latex rubber dam



Fig. 5: Preoperative clinical situation with a roeko rubber dam clamp for the anterior teeth



Fig. 6: Removal of the old restoration using a Diatech diamond FG 850L 016 12ML

ence should be given towards the use of a glass fiber reinforced or ceramic root canal post; in which a direct core build up is generally fabricated during the same appointment. Ceramic root canal posts can however also be combined with a ceramic core using the indirect technique.

There are conflicting opinions regarding the necessary properties for root canal posts and core build up materials. Some authors stipulate that root canal posts and dentin should have a similar modulus of elasticity^{1, 2, 9, 14}, whereas others claim that the rigidity of root canal posts will increase the service life expectancy of the post^{1, 15}. Neither theory is supported by adequate clinical studies. In regards to direct core build up materials, amalgam is far superior in terms of its strength and dimensional stability, although it has definite disadvantages, such as discoloration of the tooth structure due to corrosion, which rules out its use in the anterior region. Composites have a high flexural strength, while glass ceramics appear to be very suitable for fabricating a core build up in the anterior region²⁰.

Glass fiber reinforced root canal posts

According to a recently released meta-analysis⁶, prefabricated glass fiber reinforced and ceramic root canal posts failed more quickly than custom casted, metallic root canal posts. However, the failure pattern of the prefabricated glass fiber reinforced root canal posts was significantly more favourable than prefabricated or custom cast metal posts. Initial clinical data supports this supposition⁸. Although the radiopacity of glass fiber reinforced root canal posts still needs improvement, retreatment in the case of a fracture or an endodontic emergency can be completed without any problem. Unfortunately, there are still no relevant long-term studies; and a projection of the clinical behavior based on the in vitro results should be treated with caution.

Clinical case presentation

A 19 year-old patient was seeking an aesthetic improvement in the upper right

central tooth (Fig. 1). During the clinical examination a horizontal fracture line was detected on the labial aspect of the tooth, which ran approx. 4 mm coronally to the gingival margin (Fig. 1). In addition, the mesiodistal width was 1 mm less compared to tooth 21. The gingival zenith of teeth 11 and 21 were at the same level. A sufficient root canal filling on tooth 11 was visible on the radiograph (Fig. 2).

There were two alternative treatment options: insert a glass fiber reinforced root canal post and crown the tooth with an all-ceramic restoration; or perform internal bleaching and insert a glass fiber reinforced post without fitting a crown. The patient agreed to the first treatment option. The tooth shade was selected using a standardized shade guide (Fig. 3). A glass fiber reinforced root canal post was then fitted, adhesively cemented and a direct core build up fabricated using the ParaBond/ParaCore System (Fig. 4 – 28). A dental radiograph was taken afterwards to check the post (Fig. 29).

The tooth was then prepared. The sulcus was widened using the double cord



Fig. 7: Tooth with the old restoration removed and exposure of the root canal filling



Fig. 8: Removal of the root canal filling at the entrance of the root canal using a Diatech diamond FG 850L 016 12ML



Fig. 9: Removal of the root canal filling and preparation of the root canal using a Gates 1 + 2



Fig. 10: Removal of the root canal filling and preparation of the root canal using a Moser 1-3



Fig. 11: The tooth is then prepared using ParaPost drills in sequentially larger sizes until the predetermined diameter and depth is achieved.



Fig. 12: Trial placement of the ParaPost Fiber Lux

technique. The double mix technique using an addition-cured silicone was used for taking the impression (Fig. 30 – 37). In this case, AFFINIS PRECIOUS was selected, which features optimal surface affinity. This property ensures that the correction material covers the tooth surfaces immediately, even in a moist environment; and is therefore crucial for producing accurate impressions. Silver pigmentation allows excellent detail readability for assessing the quality of the impression. A direct temporary restoration was then fabricated and fitted using a silicone-based temporary cement (Fig. 38 + 39).

10 days later, the condition of the soft tissue was excellent (Fig. 40–43). The emergence profile of the completed glass ce-

ramic crown was very successful and corresponds well with the adjacent tooth (Fig. 44).

Following permanent adhesive cementation using the ParaBond/ParaCore System, the glass ceramic crown had a very acceptable length-width ratio; and the surface texture as well as the reflexion lines were an excellent match to the adjacent tooth. The gingival zenith and formation of the central papilla were highly satisfactory (Fig. 45 – 59).

Seven days after the crown was permanently fitted, the adjacent teeth were rehydrated again and exhibited a harmonious shade match with the restored tooth (Fig. 60). An excellent marginal seal

was confirmed on a follow up radiograph (Fig. 61).

Conclusion

The Monoblock Technique using ParaBond and ParaCore saves time and material; and is very versatile in an important area of restorative dentistry.

CONTACT

Priv. Doz. Dr. med. dent. Stefan J. Paul
Stadelhoferstrasse 33
CH-8001 Zürich
eMail: office@drpaul.ch



Fig. 13: Seated ParaPost Fiber Lux



Fig. 14: Irrigation of the root canal with sodium hypochlorite



Fig. 15: Drying of the root canal with sterile paper points



Fig. 16: Non-Rinse conditioner is massaged into the root canal and onto the contact surface for 30 seconds



Fig. 17: Excess Non-Rinse conditioner is removed from the root canal using a sterile paper point



Fig. 18: The tooth is then dried for 2 seconds using a gentle air stream



Fig. 19: The adhesive A+B, mixed to a 1:1 ratio, is left in the root canal and on the contact surface for 30 seconds



Fig. 20: Removal of excess adhesive from the root canal using a sterile paper point



Fig. 21: The tooth is then dried for 2 seconds using a gentle air stream



Fig. 22: Application of ParaCore core & resin cement directly into the root canal using the root canal tip



Fig. 23: Untreated ParaPost Fiber Lux Post is pre-coated with ParaCore and cemented into the root canal



Fig. 24: Removal of excess ParaCore



Fig. 25: Translucent ParaPost Fiber Lux Post is light cured for 20 s using the Coltolux LED to fixate it into place



Fig. 26: Free-hand core build up using ParaCore core & resin cement



Fig. 27: The core build up is then contoured manually.



Fig. 28: Each side of the core build up is polymerized for 20 seconds



Fig. 29: Radiograph after cementation of the post



Fig. 30: A Comprecord retraction cord size 0 is placed



Fig. 31: Preparation of the tooth using different types of Diatech diamonds



Fig. 32: Completed tooth preparation



Fig. 33: Closed gingival sulcus



Fig. 34: Second retraction cord for gingival compression



Fig. 35: Removal of the second retraction cord before taking the impression



Fig. 36: Open gingival sulcus




Fig. 37: Double mix impression using AFFINIS heavy body and AFFINIS PRECIOUS light body



Fig. 38: Trial placement of the temporary restoration fabricated using CoolTemp Natural



Fig. 39: Temporary restoration is cemented using TempoSIL 2



Fig. 40: Removal of the temporary restoration at the second appointment



Fig. 41: Prepared tooth and healed gingiva



Fig. 42: The prepared tooth is cleaned using a fluoride-free cleaning paste



Fig. 43: Prepared tooth before placement of the permanent restoration



Fig. 44: Trial placement of the permanent restoration



Fig. 45: The gingiva is slightly compressed using a retraction cord to ensure for optimal cementation



Fig. 46: Comprecord retraction cord in the sulcus



Fig. 47: The restoration is tried in again with the retraction cord in place to ensure for an accurate fit




Fig. 48: The inside of the restoration is etched...




Fig. 49: ... and silanized – always according to the manufacturer's instructions



Fig. 50: The non-rinse conditioner is massaged in for 30 seconds




Fig. 51: The non-rinse conditioner is dried using a gentle stream of air



Fig. 52: The mixed adhesive is applied onto the prepared tooth and left for 30 seconds




Fig. 53: The adhesive is dried for 2 seconds using a gentle stream of air



Fig. 54: The Root Canal Tip can be shortened using a scalpel for easy extrusion



Fig. 55: ParaCore is applied directly into the crown



Fig. 56: Initial removal of the excess cement using a sponge pellet



Fig. 57: Removal of excess cement interproximally using dental floss




Fig. 58: ParaCore can be chemically cured or light cured



Fig. 59: Occlusion is checked using Hanel articulating paper



Fig. 60: Postoperative clinical situation



Fig. 61: Postoperative radiograph

References

- 1 Anusavice KJ. Dental ceramics and metal ceramics. In: Okabe T, Takahashi S (eds). Transactions of the International Congress on Dental Ceramics, South Carolina, 1989. Academy of Dental Materials, 1989:159-172.
- 2 Assif D, Oren E, Marshak BL, Aviv I. Photoelastic analysis of stress transfer by endodontically treated teeth to the supporting structure using different restorative techniques. *J Prosthet Dent* 1989;61:535-543.
- 3 Caputo AA, Standlee JP. Restoration of endodontically treated teeth. In: Caputo AA, Standlee JP. Biomechanics in clinical dentistry. Chicago: Quintessence, 1987:185-203.
- 4 Christensen GJ. Posts, cores and patient care. *J Am Dent Assoc* 1993;124:86-90.

- ⁵ Fernandes AS, Dessai GS. Factors affecting the fracture resistance of post-core reconstructed teeth: A review. *Int J Prosthodont* 2001;14:355-363.
- ⁶ Fokkinga WA, Kreulen CM, Vallittu PK, Creugers NH. A structured analysis of in vitro failure loads and failure modes of fiber, metal, and ceramic post-and-core systems. *Int J Prosthodont* 2004;17:476-482.
- ⁷ Goracci C, Sadek FT, Fabianelli A, Tay FR, Ferrari M. Evaluation of the adhesion of fiber posts to intraradicular dentin. *Oper Dent* 2005;30:627-635.
- ⁸ Grandini S, Goracci C, Tay FR, Grandini R, Ferrari M. Clinical evaluation of the use of fiber posts and direct resin restorations for endodontically treated teeth. *Int J Prosthodont* 2005;18:399-404.
- ⁹ Hornbrook DS, Hastings JH. Use of a bondable reinforcement fiber for post and core build-up in an endodontically treated tooth: Maximizing strength and aesthetics. *Pract Periodontics Aesthet Dent* 1995;7:33-42.
- ¹⁰ Nathanson D, Ashayeri N. New aspects of restoring the endodontically treated tooth. *Alpha Omegan* 1990;83:76-80.
- ¹¹ Naumann M, Preuss A, Rosentritt M. Effect of incomplete crown ferrules on load capacity of endodontically treated maxillary incisors restored with fiber posts, composite build-ups, and all-ceramic crowns: An in vitro evaluation after chewing simulation. *Acta Odontol Scand* 2006;64:31-36.
- ¹² Paul SJ, Werder P. Clinical success of zirconium oxide posts with resin composite or glass-ceramic cores in endodontically treated teeth: A 4-year retrospective study. *Int J Prosthodont* 2004;17:524-528.
- ¹³ Randow K, Glantz PO. On cantilever loading of vital and non-vital teeth. An experimental clinical study. *Acta Odontol Scand* 1986;44:271-277.
- ¹⁴ Rudo DN, Karbhari VM. Physical behaviors of fiber reinforcement as applied to tooth stabilization. *Dent Clin North Am* 1999;43:7-35.
- ¹⁵ Scherrer SS, de Rijk WG. The fracture resistance of all-ceramic crowns on supporting structures with different elastic moduli. *Int J Prosthodont* 1993;6:462-467.
- ¹⁶ Sorensen JA, Ahn SG, Berge HX, Edelhoff D. Selection criteria for post and core materials in the restoration of endodontically treated teeth. In: *Transactions of the Conference on Scientific Criteria for Selecting Materials and Techniques in Clinical Dentistry*, Siena, Italy, September 2001. Academy of Dental Materials, 2001:67-84.
- ¹⁷ Sorensen JA, Martinoff JT. Intracoronary reinforcement and coronal coverage: A study of endodontically treated teeth. *J Prosthet Dent* 1984;51:780-784.
- ¹⁸ Stricker EJ, Göhring TN. Influence of different posts and cores on marginal adaptation, fracture resistance, and fracture mode of composite resin crowns on human mandibular premolars. An in vitro study. *J Dent* 2006;34:326-335.
- ¹⁹ Torbjörner A. Treatment management. Posts and cores. In: Karlsson S, Nilner K, Dahl BL (eds). *A textbook of fixed prosthodontics*. Stockholm: Go thia, 2000:173-186.
- ²⁰ Millar B, Sanjukta D. A comparison of marginal leakage in vitro for all-ceramic crowns luted with seven cements. King's College London Dental Institute at Guy's, King's and St. Thomas' Hospitals, London, UK. Poster Presentation, 4th ConsEuro Meeting Seville, Spain, 12. – 14. March 2009.