

REINFORCED CAD/CAM COMPOSITE FOR PERMANENT RESTORATION

BRILLIANT Crios

High performance – made brilliant



Ceramics – a smart solution
in the kitchen.

Composite – a BRILLIANT solution
in the practice.



BRILLIANT Crios

Reinforced CAD/CAM composite for permanent restoration

Composite has long been the gold standard in direct filling therapy. With BRILLIANT Crios, the well-proven dental material entered the indirect CAD/CAM restoration technology. Compared to ceramics, CAD/CAM composites offer considerable advantages.

BRILLIANT CRIOS FOR PERMANENT AESTHETIC RESTORATIONS

REINFORCED COMPOSITE

- High flexural strength for resistant restorations
- Tooth-like modulus of elasticity for a shock-absorbing effect

TOOTH-LIKE BEHAVIOR

- Blends in well to provide natural aesthetics
- High wear resistance and low abrasion on antagonist

EFFICIENT HANDLING

- No firing process required
- Can be modified and repaired
- Highly accurate grinding precision
- Effortless polishing for fast gloss

RELIABLE LUTING SYSTEM

- Secure bonding due to ONE COAT 7 UNIVERSAL
- Compatible luting material for all stated indications

SHAPES

Block	12	12 × 10 × 15 (mm)
	14	14 × 12 × 18 (mm)
Disc	H14	14 × 98.5 (mm)
	H18	18 × 98.5 (mm)

SHADES

Low Translucent (LT)

Bleach, A1, A2, A3, A3.5, B1, B2, B3, C2

High Translucent (HT)

A1, A2, A3, B1

Super Translucent (ST)

Bleach, Universal


NEW

INDICATIONS


BRILLIANT Crios is the ideal choice for single tooth restorations in the anterior as well as posterior region of teeth. This includes all conventional indications such as **inlays, onlays, crowns and veneers**.



Inlay




Onlay



Crown



Veneer



REINFORCED COMPOSITE

The outstanding mechanical properties of BRILLIANT Crios are the result of controlled, stress-free thermal curing. The multimodal composition of dental glass and amorphous silica in combination with a reinforcing resin matrix make BRILLIANT Crios the ideal material for permanent single tooth restorations.

A STRONG MATERIAL

The quality of BRILLIANT Crios is also reflected by its high flexural strength. The compact structure increases the material resistance and reduces the risk of fracturing.

Source: internal data | * not a trademark of COLTENE


SHOCK-ABSORBING EFFECT

The dentin-like modulus of elasticity of BRILLIANT Crios reduces stress peaks during chewing, resulting in decreased risk of material fatigue. Due to its shock-absorbing effect, it is well suited for implant restorations. As tooth substance and BRILLIANT Crios have a similar modulus of elasticity, loading stress is considerably reduced at the margins.

Source: R. Böhner et al., #597, IADR 2015, Antalya, Turkey
* not a trademark of COLTENE


BIAXIAL FLEXURAL STRENGTH

Measured in MPa




MODULUS OF ELASTICITY

Measured in GPa



TOOTH-LIKE BEHAVIOR

TWO-BODY WEAR



Source: Stawarczyk et al., J Mech Behav Biomed Mater 55, 1-11 (2015) | * not a trademark of COLTENE

WEAR RESISTANT AND GENTLE AT THE SAME TIME


Studies confirm the excellent abrasion resistance of BRILLIANT Crios. Nonetheless, the wear-resistant composite material proves to be especially kind to the antagonist teeth in comparison to ceramic restorations. BRILLIANT Crios is suitable as a permanent restoration material – even under considerable chewing forces – and thus offers the best of both worlds.



NATURAL AESTHETICS WITH A GLOSSY RESULT

BRILLIANT Crios blends in well and conveys natural aesthetics. A broad spectrum of 15 shades in three translucencies perfectly covers the needs in the anterior as well as posterior region of teeth. Immediately after grinding, BRILLIANT Crios displays a homogeneous surface structure. The restoration can therefore be brought to a high gloss very quickly. Easy repolishing, for example during prophylaxis, is possible at any time if required.


EFFICIENT HANDLING




EXCEPTIONAL GRINDING PRECISION

BRILLIANT Crios is distinguished by superior precision. Even when ceramic materials have already start chipping, BRILLIANT Crios can still be ground with ease and precision in the case of extremely thin margins. This exceptional grinding precision offers greater freedom of preparation.


BRILLIANT Crios




Lava Ultimate*



VITA Enamic*



IPS Empress CAD*



Ground tip 0.1 mm (Sirona inLab MC XL)

Source: C. Kopfmann et al., Eur. J. Prosthodont. Rest. Dent., EMDC Special, P21, 2015, Nuremberg | * not a trademark of COLTENE

TIME AND COST EFFICIENT

- Effortless polishing for fast gloss
- No firing process required
- Quick grinding at high precision
- Easy to modify and repair

Photo: Dr. Paulo Monteiro



CLINICAL CASE

Restored teeth

26 with BRILLIANT Crios A2 HT

27 with BRILLIANT Crios A2 LT

Luting

ONE COAT 7 UNIVERSAL

with BRILLIANT EverGlow A2/B2

1. Initial situation



2. Prepare



3. CAD/CAM restoration (27)



4. Sandblasting



5. Bonding




6. Etching



7. Bonding




8. Light-curing



9. Luting



10. Light-curing




11. Final result



RELIABLE LUTING SYSTEM

As part of its adhesive filling therapy, COLTENE offers a complementary proven luting system for BRILLIANT Crios. The key component of the luting concept is the universal adhesive ONE COAT 7 UNIVERSAL. The smooth BRILLIANT EverGlow composite or the dual-curing resin cement SoloCem offer

indication-appropriate solutions, depending on the wall thickness of the restoration and the luting substrate, ranging from natural dentin to a metal abutment.



BRILLIANT EverGlow

Light-curing composite

Ideal for luting indirect restorations with thin wall design to natural tooth substance or composite:

- Smooth consistency and excellent wettability to simplify restoration placement
- Exceptional aesthetic and long-lasting results

For restorations up to 3 mm wall thickness

SoloCem

Self-adhesive, dual-curing resin cement

Convincing universal luting solution for any kind of substrate:

- To tooth/composite using ONE COAT 7 UNIVERSAL as additional bonding agent
- To metal/ceramic abutments, without requiring an additional bonding step

For restorations up to 5 mm wall thickness

SECURE BONDING

ONE COAT 7 UNIVERSAL


ONE COAT 7 UNIVERSAL is matched perfectly to the BRILLIANT Crios components, resulting in strong anchorage and a reliable adhesion. This is why the ONE COAT 7 UNIVERSAL adhesive is to be used in conjunction with BRILLIANT Crios.

 BRILLIANT Crios
 ONE COAT 7 UNIVERSAL

The monomers (red) contained in ONE COAT 7 UNIVERSAL penetrate into the polymerized resin matrix of BRILLIANT Crios




Polymer chains are formed during polymerization of the monomers and these are anchored in the resin matrix.



The double bonds contained in the resin matrix form chains during polymerization of the monomers.

SHEAR BOND STRENGTH WITH ONE COAT 7 UNIVERSAL

Measured in MPa



The BRILLIANT EverGlow and SoloCem luting materials demonstrate excellent bonding strengths to give a secure adhesion to BRILLIANT Crios.

Source: Internal data
According to Watanabe

FURTHER PROPERTIES

ONE COAT 7 UNIVERSAL is a light-curing one component universal adhesive.

- Exceptional bonding values on enamel and dentin
- Excellent marginal seal
- Minimized postoperative sensitivity
- For all application techniques (Self-Etch, Selective Etch, Total Etch)



ORDER INFORMATION

BRILLIANT Crios Kits

60020095 Intro Kit CEREC 12

Each 1 × A2 LT 12, A3 LT 12, A3.5 LT 12, A2 HT 12, A3 HT 12
5 × 0.1 ml ONE COAT 7 UNIVERSAL Single Dose

60020097 Intro Kit CEREC 14

Each 1 × A2 LT 14, A3 LT 14, A3.5 LT 14, A2 HT 14, A3 HT 14
5 × 0.1 ml ONE COAT 7 UNIVERSAL Single Dose

60020099 Intro Kit Planmill 14

Each 1 × A2 LT 14, A3 LT 14, A3.5 LT 14, A2 HT 14, A3 HT 14
5 × 0.1 ml ONE COAT 7 UNIVERSAL Single Dose



BRILLIANT Crios Refills

Translucency	CEREC/inLab				
	Low (LT)		High (HT)		Super (ST)
Size	12	14	12	14	14
Content	5 pcs	5 pcs	5 pcs	5 pcs	5 pcs
Bleach (BL)	60019998	60019985			60024174
A1	60019999	60019986	60020007	60019994	
A2	60020000	60019987	60020008	60019995	
A3	60020001	60019988	60020009	60019996	
A3.5	60020002	60019989			
B1	60020003	60019990	60020010	60019997	
B2	60020004	60019991			
B3	60020005	60019992			
C2	60020006	60019993			
Universal (UN)					60024175

Planmill	
Low (LT)	High (HT)
14	14
5 pcs	5 pcs
60020011	
60020012	60020020
60020013	60020021
60020014	60020022
60020015	
60020016	60020023
60020017	
60020018	
60020019	

Translucency	Disc					
	Low (LT)		High (HT)		Super (ST)	
Height (mm)	14	18	14	18	14	18
Ø (mm)	98.5	98.5	98.5	98.5	98.5	98.5
Content	1 pc	1 pc	1 pc	1 pc	1 pc	1 pc
Bleach (BL)					60024179	60024177
A1	60022899	60022886	60022907	60022894		
A2	60022900	60022887	60022908	60022895		
A3	60022901	60022888	60022909	60022896		
A3.5	60022902	60022889				
Universal (UN)					60024180	60024178



© COLTENE – www.coltene.com

Coltene/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
T +1 330 916 8800
F +1 330 645 8704
info.us@coltene.com

Coltene/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info.ch@coltene.com

Coltene/Whaledent GmbH + Co. KG **EC REP**
Raiffeisenstrasse 30
89129 Langenau / Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info.de@coltene.com

PN:31548B

COLTENE