

## BRILLIANT EverGlow Farbkombinationen

Zahnfarbe (VITA)	Universal							Translucent		Opaque				
	BL	A1/B1	A2/B2	A3/D3	A3.5/B3	C2/C3	A4/C4	BL Trans	Trans	OBL	OA1	OA2	OA3	OA4
Bleach	xx							xx		xx				
A1		xx						xx	xx		xx	x		
A2			xx					x	xx		x	xx	x	
A3				xx				x	xx			x	xx	
A3.5					xx				xx				x	x
A4							xx		xx					xx
B1		xx						xx	xx	x	x			
B2			xx					x	xx		x	x		
B3					xx			x	xx			x	x	
D3				xx					xx				x	
C2						xx			xx				x	
C3						xx			xx				x	
C4							xx		xx					x
Jung / Aufgehell														
			Erwachsene											
				Ältere										

xx besonders geeignet  
x geeignet