


POLISHING INSTRUMENTS
FOR THE DENTAL LABORATORY (HP AND UM)


RODIA TECH



ISO 658 900 303 524 ...



ISO 658 900 372 524 ...



ISO 658 900 303 514 ...



ISO 658 900 372 514 ...



CERAPRESHINE		
ISO Ø 1/10 mm	220	220
Head L in mm	4	3
Polisher type	3236UM	3237UM
UM 10 pcs	230134AA	230136AA

CERASHINE	
220	220
4	3
3336UM	3337UM
230159AA	230161AA

Ceramic ① ②

ISO 805 104 173 534 ...



ISO 805 104 107 534 ...



ISO 805 104 042 534 ...



ISO 803 104 243 524 ...



ISO 803 104 372 524 ...



ISO 803 104 303 524 ...



ISO 803 104 243 514 ...



ISO 803 104 372 514 ...



ISO 803 104 303 514 ...



ZIRCOGRIND PLUS			
ISO Ø 1/10 mm	040	050	150
Head L in mm	10	13	3
Polisher type	9640HP	9670HP	9636HP
HP 5 pcs	230373AA	230375AA	230371AA
HP 1 pc			

ZIRCOPOL PLUS		
040	110	260
13	2	2
9703HP	9731HP	9738HP
230377AA	230383AA	
		230387AA

ZIRCOSHINE PLUS		
040	110	260
13	2	2
9803HP	9831HP	9838HP
230389AA	230395AA	
		230399AA

CAD/CAM Restoration, Zirconia ① ② ③

ISO 803 104 243 534 ...



ISO 803 104 372 534 ...



ISO 803 104 303 534 ...



ISO 803 104 243 524 ...



ISO 803 104 372 524 ...



ISO 803 104 303 524 ...



ISO 803 104 372 514 ...



ISO 803 104 303 514 ...



CERAFIN PLUS			
ISO Ø 1/10 mm	040	110	260
Head L in mm	13	2	2
Polisher type	9303HP	9331HP	9338HP
HP 5 pcs	230315AA	230323AA	
HP 1 pc			230329AA

CERAPRESHINE PLUS		
040	110	260
13	2	2
9403HP	9431HP	9438HP
230333AA	230341AA	
		230347AA

CERASHINE PLUS	
110	260
2	2
9531HP	9538HP
230359AA	
	230365AA

Ceramic, Metal ① ② ③

ISO 658 900 030 533 ...



ISO 658 900 114 533 ...



ISO 658 900 372 533 ...



ISO 658 900 303 533 ...



ISO 658 104 243 523 ...



UNIPOL					UNISHINE
ISO Ø 1/10 mm	110	070	220	220	055
Head L in mm	6	20	3	4	15
Polisher type	4116UM	4171UM	4137UM	4138UM	4203HP
HP 10 pcs					230205AA
UM 10 pcs	230180AA	60019535	230196AA	60019534	
Universal					

ISO 658 104 292 513 ...



ISO 658 900 372 513 ...



ISO 658 104 292 503 ...



ISO 658 900 372 503 ...



ISO 803 104 292 523 ...



ISO 803 104 304 523 ...



PREBROWN			HIGREEN		
ISO Ø 1/10 mm	050	220	050	220	
Head L in mm	16	3	16	3	
Polisher type	1104HP	1137UM	1204HP	1237UM	
HP 10 pcs	230011AA		230038AA		
UM 10 pcs		230025AA		230052AA	
Precious Metal Alloy, Amalgam					

DIASHINE COMPOMANT PLUS		
ISO Ø 1/10 mm	050	145
Head L in mm	16	2.5
Polisher type	9104HP	9121HP
HP 5 pcs	60021794	60021795
Composite, Compomere		

ISO 653 900 372 534 ...



ISO 653 900 303 534 ...



ISO 653 900 114 534 ...



STEELMASTER			
ISO Ø 1/10 mm	220	220	060
Head L in mm	3	4	22
Polisher type	8137UM	8138UM	8171UM
UM 50 pcs	232010AA	232011AA	232012AA
Hard alloy			

ISO 658 104 274 533 ...



ISO 658 104 273 533 ...



ISO 658 104 274 523 ...



ISO 658 104 273 523 ...



ISO 658 104 274 513 ...




ACRYPREPOL, ACRYPOL, ACRYSHINE					
ISO Ø 1/10 mm	100	150	100	150	100
Head L in mm	24	17	24	17	24
Polisher type	6150HP	6160HP	6250HP	6260HP	6350HP
HP 5 pcs	230247AA	230249AA	230253AA	230255AA	230259AA
Acrylic					

Order information

ISO number
(material, shank, shape, grit size,
variable diameter)

Head diameter in 1/10 mm
Head length
Polisher description
Packaging sizes available

ISO 803 104 243 514 ...



Illustration

ISO Ø 1/10 mm	040
Head L in mm	13
Polisher type	9803HP
HP 5 pcs	230389AA
HP 1 pc	

Order numbers in respective
packaging size

Head Shapes

ISO Head Shape	Description
030	Cup
042	Wheel
107	Cylinder
114	Cylinder
173	Taper
243	Point
273	Dome
274	Bullet
292	Torpedo
303	Large Lens
304	Small Lens
372	Disc

Recommended Speed Ranges

Polishing System	RPM
CERAPRESHINE CERASHINE	5'000–10'000
ZIRCOGRIND PLUS	8'000–12'000
ZIRCOPOL PLUS 9738HP	7'000–10'000 5'000–7'000
ZIRCOSHINE PLUS 9838HP	7'000–10'000 5'000–7'000
CERAFIN PLUS 9338HP	7'000–10'000 5'000–7'000
CERAPRESHINE PLUS 9438HP	7'000–10'000 5'000–7'000
CERASHINE PLUS 9538HP	7'000–10'000 5'000–7'000

Polishing System	RPM
UNIPOL UNISHINE	5'000–10'000
PREBROWN HIGREEN	5'000–10'000
STEELMASTER	5'000–10'000
ACRYPREPOL ACRYPOL ACRYSHINE	5'000–7'000
DIASHINE COMPOMANT PLUS	7'000–10'000 Prepolishing 3'000–8'000 High lustre

CH Coltene/Whaledent AG
Feldwiesenstr. 20
9450 Altstätten / Switzerland
T + 41 71 757 53 00
F + 41 71 757 53 01
info.ch@coltene.com

DE Coltene/Whaledent GmbH + Co. KG
Raiffeisenstraße 30
89129 Langenau / Germany
T + 49 7345 805 0
F + 49 7345 805 201
info.de@coltene.com

US Coltene/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
T + 1 330 916 8800
F + 1 330 916 7077
info.us@coltene.com

© 2016 Coltene/Whaledent AG – www.coltene.com



COLTENE

