

Alpen Conversion List

Important information **EN**

We are delighted to inform our valued clients that COLTENE has made significant improvements regarding the packaging of our Alpen Diamonds:

- User-friendly blister packaging format
- Securely sealed, preserving integrity of the product
- Equipped with an easy peel feature

Additionally, we would like to notify you about certain changes in ring colors and grit sizes to align with the ISO 7711 standard (see table below):

- The black “super coarse” (MLX) has been reclassified as green “coarse” (C)
- The green “coarse” (ML) has been reclassified as blue “medium” (M)

Furthermore, there have been modifications to the letters in the names of our burs:

- The black MLX instruments have been transformed to green C
- The blue ML to blue M
- The yellow XF has been replaced by yellow FF



Please note that there have been changes to the item numbers. For your reference, we have included a conversion table that lists the old and new item numbers. This conversion list only includes new items, that have an equivalent old item. However, the new portfolio also includes completely new items, that are not shown in this list. You can access the PDF version of this table through the provided link: <https://pxl.coltene.com/alpenqr> or by scanning the QR code.

We appreciate your understanding as we strive to ensure precision and quality and provide the best possible solutions for you.

Your COLTENE Team

| Old color ring - index (grit type) | Old mean grit value | New color ring - index (grit type) | New grit value range and mean |
|------------------------------------|---------------------|---|---------------------------------|
| Black - SC (super coarse) | 140 µm | Green - C (coarse) | 107-181 µm; Mean 151 µm |
| Green - C (coarse) | 120 µm | Blue - M (medium) | 64-126 µm; Mean 107 µm |
| Grey - M (medium) | 105 µm | Blue - M (medium) | 64-126 µm; Mean 107 µm |
| Red - F (fine) | 40 µm | Red - F (fine) | 27-76 µm; Mean 46 µm |
| Yellow - EF (extra fine) | 30 µm | Yellow - EF (extra fine) | 10-36 µm; Mean 25 µm |
| White - UF (ultra fine) | 20 µm | not available anymore | 4-14 µm; Mean 8 µm |
| Black - TSX (Turbo Cut) | 140µm | Green - C (coarse) or Blue - M (medium) | does not exist according to ISO |

| REF new | Description new | REF old | Description old |
|----------|--------------------------------------|---------------|---------------------|
| 60031753 | G368-314-016-03.5-EF, Alpen Diamonds | R368F016FG | G368-314-016-3.5-XF |
| 60031755 | G368-314-016-03.5-F, Alpen Diamonds | R368F016FG | G368-314-016-3.5-F |
| 60031757 | G368-314-016-03.5-M, Alpen Diamonds | R368C016FG | G368-314-016-3.5-ML |
| 60031759 | G368-314-020-05.0-C, Alpen Diamonds | R368SC020FG | G368-314-020-5-MLX |
| 60031761 | G368-314-020-05.0-EF, Alpen Diamonds | R368F020FG | G368-314-020-5-F |
| 60031763 | G368-314-020-05.0-F, Alpen Diamonds | R368F020FG | G368-314-020-5-F |
| 60031767 | G368-314-020-05.0-M, Alpen Diamonds | R368C020FG | G368-314-020-5-ML |
| 60031769 | G368-314-023-05.0-C, Alpen Diamonds | R368SC023FG | G368-314-023-5-MLX |
| 60031771 | G368-314-023-05.0-EF, Alpen Diamonds | R368F023FG | G368-314-023-5-F |
| 60031773 | G368-314-023-05.0-F, Alpen Diamonds | R368F023FG | G368-314-023-5-F |
| 60031777 | G368-314-023-05.0-M, Alpen Diamonds | R368C023FG | G368-314-023-5-ML |
| 60031811 | G379-314-016-03.5-EF, Alpen Diamonds | R379F016FG | G379-314-016-3.5-XF |
| 60031783 | G379-314-016-03.5-F, Alpen Diamonds | R379F016FG | G379-314-016-3.5-F |
| 60031785 | G379-314-016-03.5-M, Alpen Diamonds | R379C016FG | G379-314-016-3.5-ML |
| 60031787 | G379-314-018-03.5-EF, Alpen Diamonds | R379F018FG | G379-314-018-3.5-XF |
| 60031789 | G379-314-018-03.5-F, Alpen Diamonds | R379F018FG | G379-314-018-3.5-F |
| 60031793 | G379-314-018-03.5-M, Alpen Diamonds | R379C018FG | G379-314-018-3.5-ML |
| 60031795 | G379-314-023-05.0-C, Alpen Diamonds | R379SC023FG | G379-314-023-5-MLX |
| 60031797 | G379-314-023-05.0-EF, Alpen Diamonds | R379F023FG | G379-314-023-5-F |
| 60031801 | G379-314-023-05.0-F, Alpen Diamonds | R379F023FG | G379-314-023-5-F |
| 60031807 | G379-314-023-05.0-M, Alpen Diamonds | R379C023FG | G379-314-023-5-ML |
| 60031811 | G379T-314-023-05.0-C, Alpen Diamonds | R379T023FG | G379-314-023-5-TSX |
| 60031813 | G379Z-314-023-05.0-M, Alpen Diamonds | ZREX379M023FG | ZR379-314-023-5-M |
| 60031815 | G390-314-014-03.0-EF, Alpen Diamonds | R390F014FG | G390-314-014-3-F |
| 60031817 | G390-314-014-03.0-F, Alpen Diamonds | R390F014FG | G390-314-014-3-F |
| 60031819 | G392-314-016-06.0-EF, Alpen Diamonds | R392F016FG | G392-314-016-6-F |
| 60031821 | G392-314-016-06.0-F, Alpen Diamonds | R392F016FG | G392-314-016-6-F |
| 60031825 | G392-314-016-06.0-M, Alpen Diamonds | R392C016FG | G392-314-016-6-ML |

| REF new | Description new | REF old | Description old |
|----------|--------------------------------------|---------------|-----------------------|
| 60031827 | G392-314-016-08.0-F, Alpen Diamonds | R392F016FG | G392-314-016-8-F |
| 60031829 | G392-314-016-08.0-M, Alpen Diamonds | R392C016FG | G392-314-016-8-ML |
| 60031831 | G392-314-018-08.0-EF, Alpen Diamonds | R392F018FG | G392-314-018-8-F |
| 60031833 | G392-314-018-08.0-F, Alpen Diamonds | R392F018FG | G392-314-018-8-F |
| 60031835 | G392-314-018-08.0-M, Alpen Diamonds | R392C018FG | G392-314-018-8-ML |
| 60031837 | G392-314-021-10.0-F, Alpen Diamonds | R392F021FG | G392-314-021-10-F |
| 60031839 | G392-314-021-10.0-M, Alpen Diamonds | R392C021FG | G392-314-021-10-ML |
| 60031841 | G6051-314-018-06.0-M, Alpen Diamonds | R6051C018FG | G6051-314-018-6-ML |
| 60031843 | G6052-314-018-08.0-C, Alpen Diamonds | R6051C018FG | G6051-314-018-6-MLX |
| 60031845 | G6052-314-018-08.0-M, Alpen Diamonds | R6052C018FG | G6052-314-018-8-ML |
| 60031847 | G6055-314-018-08.0-C, Alpen Diamonds | R6055C018FG | G6055-314-018-8-MLX |
| 60031849 | G801-314-009-00.9-M, Alpen Diamonds | R801C009FG | G801-314-009-ML |
| 60031851 | G801-314-010-01.0-M, Alpen Diamonds | R801C010FG | G801-314-010-ML |
| 60031853 | G801-314-012-01.2-M, Alpen Diamonds | R801C012FG | G801-314-012-ML |
| 60031857 | G801-314-014-01.4-F, Alpen Diamonds | R801F014FG | G801-314-014-F |
| 60031859 | G801-314-014-01.4-M, Alpen Diamonds | R801C014FG | G801-314-014-ML |
| 60031865 | G801-314-016-01.6-M, Alpen Diamonds | R801C016FG | G801-314-016-ML |
| 60031871 | G801-314-018-01.8-F, Alpen Diamonds | R801F018FG | G801-314-018-F |
| 60031873 | G801-314-018-01.8-M, Alpen Diamonds | R801C018FG | G801-314-018-ML |
| 60031877 | G801-314-020-02.0-F, Alpen Diamonds | R801F021FG | G801-314-020-F |
| 60031879 | G801-314-020-02.0-M, Alpen Diamonds | R801C021FG | G801-314-020-ML |
| 60031881 | G801-314-023-02.3-C, Alpen Diamonds | R801C023FG | G801-314-023-MLX |
| 60031883 | G801-314-023-02.3-EF, Alpen Diamonds | R801F023FG | G801-314-023-XF |
| 60031885 | G801-314-023-02.3-F, Alpen Diamonds | R801C023FG | G801-314-023-F |
| 60031887 | G801-314-023-02.3-M, Alpen Diamonds | R801C023FG | G801-314-023-ML |
| 60031891 | G801-314-027-02.7-M, Alpen Diamonds | R801C027FG | G801-314-027-ML |
| 60031893 | G801-314-029-02.9-F, Alpen Diamonds | R801F029FG | G801-314-029-F |
| 60031895 | G801-314-029-02.9-M, Alpen Diamonds | R801C029FG | G801-314-029-ML |
| 60031897 | G801M-314-009-00.9-F, Alpen Diamonds | R801MF009FG | G801M-314-009-0.6-F |
| 60031899 | G801Z-314-012-01.2-M, Alpen Diamonds | ZREX801M012FG | ZR801-314-012-M |
| 60031901 | G801Z-314-014-01.4-M, Alpen Diamonds | ZREX801M014FG | ZR801-314-014-M |
| 60031903 | G801Z-314-018-01.8-M, Alpen Diamonds | ZREX801M018FG | ZR801-314-018-M |
| 60031911 | G805-314-012-00.9-M, Alpen Diamonds | R805C012FG | G805-314-012-0.9-ML |
| 60031913 | G805-314-014-01.4-M, Alpen Diamonds | R805C014FG | G805-314-014-1.4-ML |
| 60031915 | G805-314-016-01.4-M, Alpen Diamonds | R805C016FG | G805-314-016-1.4-ML |
| 60031919 | G806-314-012-02.2-M, Alpen Diamonds | R806C012FG | G806-314-012-2.2-ML |
| 60031921 | G806-314-014-02.2-M, Alpen Diamonds | R806C014FG | G806-314-014-2.2-ML |
| 60031923 | G806-314-016-02.5-M, Alpen Diamonds | R806C016FG | G806-314-016-2.5-ML |
| 60031927 | G807-314-012-03.5-M, Alpen Diamonds | R807C012FG | G807-314-012-3.5-ML |
| 60031929 | G807-314-016-04.0-M, Alpen Diamonds | R807C016FG | G807-314-016-4-ML |
| 60031931 | G807-314-018-05.0-M, Alpen Diamonds | R807C018FG | G807-314-018-5-ML |
| 60031933 | G811-314-031-04.2-M, Alpen Diamonds | R811C031FG | G811-314-031-4.2-ML |
| 60031935 | G811-314-033-05.0-C, Alpen Diamonds | R811C033FG | G811-314-033-5-MLX |
| 60031937 | G811-314-033-05.0-M, Alpen Diamonds | R811C033FG | G811-314-033-5-ML |
| 60031939 | G811-314-037-07.0-C, Alpen Diamonds | R811C037FG | G811-314-037-7-MLX |
| 60031941 | G811-314-037-07.0-M, Alpen Diamonds | R811C037FG | G811-314-037-7-ML |
| 60031943 | G811-314-048-08.0-M, Alpen Diamonds | R811C048FG | G811-314-048-8-ML |
| 60031945 | G811T-314-033-05.0-C, Alpen Diamonds | R811TC033FG | G811-314-033-5-TSX |
| 60031947 | G813-314-014-01.5-M, Alpen Diamonds | R813C014FG | G813-314-014-1.5-ML |
| 60031951 | G813-314-018-02.0-M, Alpen Diamonds | R813C018FG | G813-314-018-2-ML |
| 60031955 | G818-314-050-00.7-M, Alpen Diamonds | R818C050FG | G818-314-050-0.7-ML |
| 60031961 | G828-314-030-01.0-M, Alpen Diamonds | R828C030FG | G828-314-030-1.0-ML |
| 60031963 | G830-313-010-02.7-M, Alpen Diamonds | R830C010FG | G830-313-010-2.7-MLSS |
| 60031967 | G830-314-010-02.7-M, Alpen Diamonds | R830C010FG | G830-314-010-2.7-ML |
| 60031969 | G830-314-012-02.7-M, Alpen Diamonds | R830C012FG | G830-314-012-2.7-ML |
| 60031973 | G830-314-016-02.7-M, Alpen Diamonds | R830C016FG | G830-314-016-2.7-ML |
| 60031981 | G830L-314-012-04.0-M, Alpen Diamonds | R830L012FG | G830L-314-012-4-ML |
| 60031985 | G830L-314-014-05.0-M, Alpen Diamonds | R830L014FG | G830L-314-014-5-ML |
| 60031987 | G830L-314-016-05.0-C, Alpen Diamonds | R830L016FG | G830L-314-016-5-ML |
| 60031995 | G834-314-016-06.8-M, Alpen Diamonds | R834C016FG | G834-314-016-6.8-ML |
| 60031997 | G834-314-021-06.8-M, Alpen Diamonds | R834C021FG | G834-314-021-6.8-ML |
| 60031999 | G835-313-010-04.0-M, Alpen Diamonds | R835C010SS | G835-313-010-4-MLSS |
| 60032001 | G835-314-008-03.0-M, Alpen Diamonds | R835C008FG | G835-314-008-3-ML |
| 60032003 | G835-314-010-04.0-M, Alpen Diamonds | R835C010FG | G835-314-010-4-ML |
| 60032007 | G835-314-012-04.0-M, Alpen Diamonds | R835C012FG | G835-314-012-4-ML |
| 60032009 | G835-314-014-04.0-M, Alpen Diamonds | R835C014FG | G835-314-014-4-ML |
| 60032013 | G835R-313-010-04.0-M, Alpen Diamonds | R835RC010SS | G835R-313-010-4-MLSS |
| 60032015 | G835R-314-008-03.0-M, Alpen Diamonds | R835RC008FG | G835R-314-008-3-ML |
| 60032017 | G835R-314-010-04.0-M, Alpen Diamonds | R835RC010FG | G835R-314-010-4-ML |
| 60032021 | G835R-314-012-04.0-M, Alpen Diamonds | R835RC012FG | G835R-314-012-4-ML |
| 60032023 | G835R-314-014-04.0-M, Alpen Diamonds | R835RC014FG | G835R-314-014-4-ML |
| 60032025 | G836-314-012-06.0-F, Alpen Diamonds | R836F012FG | G836-314-012-6-F |
| 60032027 | G836-314-012-06.0-M, Alpen Diamonds | R836C012FG | G836-314-012-6-ML |
| 60032029 | G836-314-014-06.0-C, Alpen Diamonds | R836SC014FG | G836-314-014-6-MLX |
| 60032031 | G836-314-014-06.0-M, Alpen Diamonds | R836C014FG | G836-314-014-6-ML |
| 60032037 | G836R-314-010-06.0-M, Alpen Diamonds | R836RC010FG | G836R-314-010-6-ML |
| 60032039 | G836R-314-012-06.0-M, Alpen Diamonds | R836RC012FG | G836R-314-012-6-ML |
| 60032041 | G836R-314-014-06.0-F, Alpen Diamonds | R836RC014FG | G836R-314-014-6-F |
| 60032043 | G836R-314-014-06.0-M, Alpen Diamonds | R836RC014FG | G836R-314-014-6-ML |

| REF new | Description new | REF old | Description old |
|----------|--------------------------------------|-------------|---------------------|
| 60032045 | G836T-314-014-06.0-M, Alpen Diamonds | R836TC014FG | G836-314-014-6-TSX |
| 60032047 | G836T-314-016-06.0-C, Alpen Diamonds | R836TC016FG | G836-314-016-6-TSX |
| 60032049 | G837-314-012-07.0-M, Alpen Diamonds | R837C012FG | G837-314-012-7-ML |
| 60032051 | G837-314-012-08.0-M, Alpen Diamonds | R837C012FG | G837-314-012-8-ML |
| 60032053 | G837-314-014-08.0-M, Alpen Diamonds | R837C014FG | G837-314-014-8-ML |
| 60032055 | G837-314-016-08.0-C, Alpen Diamonds | R837C016FG | G837-314-016-8-MLX |
| 60032057 | G837-314-016-08.0-M, Alpen Diamonds | R837C016FG | G837-314-016-8-ML |
| 60032061 | G837R-314-012-08.0-M, Alpen Diamonds | R837RC012FG | G837R-314-012-8-ML |
| 60032063 | G837R-314-014-08.0-F, Alpen Diamonds | R837RF014FG | G837R-314-014-8-F |
| 60032065 | G837R-314-014-08.0-M, Alpen Diamonds | R837RC014FG | G837R-314-014-8-ML |
| 60032067 | G837R-314-016-08.0-M, Alpen Diamonds | R837RC016FG | G837R-314-016-8-ML |
| 60032069 | G837T-314-014-08.0-M, Alpen Diamonds | R837TC014FG | G837-314-014-8-TSX |
| 60032071 | G837T-314-016-08.0-C, Alpen Diamonds | R837TC016FG | G837-314-016-8-TSX |
| 60032075 | G838-314-010-04.0-M, Alpen Diamonds | R838C010FG | G838-314-010-4-ML |
| 60032077 | G838-314-012-04.0-M, Alpen Diamonds | R838C012FG | G838-314-012-4-ML |
| 60032081 | G838M-314-007-02.5-F, Alpen Diamonds | R838MF002FG | G838M-314-007-2.4-F |
| 60032083 | G839-314-010-00.2-M, Alpen Diamonds | R839C010FG | G839-314-010-ML |
| 60032085 | G839-314-012-00.2-F, Alpen Diamonds | R839F012FG | G839-314-012-F |
| 60032087 | G839-314-012-00.2-M, Alpen Diamonds | R839C012FG | G839-314-012-ML |
| 60032089 | G839-314-014-00.2-F, Alpen Diamonds | R839F014FG | G839-314-014-F |
| 60032091 | G839-314-014-00.2-M, Alpen Diamonds | R839C014FG | G839-314-014-ML |
| 60032095 | G845-314-010-04.0-M, Alpen Diamonds | R845C010FG | G845-314-010-4-ML |
| 60032097 | G845-314-012-04.0-M, Alpen Diamonds | R845C012FG | G845-314-012-4-ML |
| 60032103 | G845-314-016-04.0-M, Alpen Diamonds | R845C016FG | G845-314-016-4-ML |
| 60032105 | G845R-314-016-04.0-M, Alpen Diamonds | R845RC016FG | G845R-314-016-4-ML |
| 60032107 | G845R-314-025-04.0-F, Alpen Diamonds | R845RF025FG | G845R-314-025-4-F |
| 60032109 | G845R-314-025-04.0-M, Alpen Diamonds | R845RC025FG | G845R-314-025-4-ML |
| 60032111 | G846-314-016-06.0-M, Alpen Diamonds | R846C016FG | G846-314-016-6-ML |
| 60032113 | G846-314-016-07.0-M, Alpen Diamonds | R846C016FG | G846-314-016-7-ML |
| 60032115 | G846-314-018-07.0-M, Alpen Diamonds | R846C018FG | G846-314-018-7-ML |
| 60032117 | G846-314-023-07.0-M, Alpen Diamonds | R846C023FG | G846-314-023-7-ML |
| 60032119 | G846R-314-016-06.0-F, Alpen Diamonds | R846RF016FG | G846R-314-016-6-F |
| 60032121 | G846R-314-016-06.0-M, Alpen Diamonds | R846RC016FG | G846 |

| REF new | Description new | REF old | Description old |
|----------|---------------------------------------|---------------|----------------------|
| 60032227 | G850L-314-016-12.0-M, Alpen Diamonds | R850LC016FG | G850L-314-016-12-ML |
| 60032229 | G850L-314-018-12.0-M, Alpen Diamonds | R850LC018FG | G850L-314-018-12-M |
| 60032233 | G850T-314-014-10.0-M, Alpen Diamonds | R850TC014FG | G850-314-014-10-TSX |
| 60032235 | G850T-314-016-10.0-C, Alpen Diamonds | R850TC016FG | G850-314-016-10-TSX |
| 60032237 | G850T-314-018-08.0-C, Alpen Diamonds | R850TC018FG | G850-314-018-8-TSX |
| 60032239 | G850T-314-018-10.0-C, Alpen Diamonds | R850TC018FG | G850-314-018-10-TSX |
| 60032241 | G850Z-314-016-10.0-M, Alpen Diamonds | ZREX850M016FG | ZREX850-314-016-10-M |
| 60032243 | G850Z-314-018-10.0-M, Alpen Diamonds | ZREX850M018FG | ZREX850-314-018-10-M |
| 60032245 | G851-314-012-08.0-M, Alpen Diamonds | R851C012FG | G851-314-012-8-ML |
| 60032247 | G851-314-016-08.0-M, Alpen Diamonds | R851C016FG | G851-314-016-8-ML |
| 60032249 | G852-314-012-06.0-EF, Alpen Diamonds | R852EF012FG | G852-314-012-6-F |
| 60032251 | G852-314-012-06.0-F, Alpen Diamonds | R852F012FG | G852-314-012-6-F |
| 60032253 | G852-314-012-06.0-M, Alpen Diamonds | R852C012FG | G852-314-012-6-ML |
| 60032255 | G852-314-016-06.0-EF, Alpen Diamonds | R852EF016FG | G852-314-016-6-F |
| 60032259 | G852-314-016-06.0-M, Alpen Diamonds | R852C016FG | G852-314-016-6-ML |
| 60032261 | G853-314-008-03.5-EF, Alpen Diamonds | R853EF008FG | G853-314-008-3.5-F |
| 60032263 | G853-314-008-03.5-M, Alpen Diamonds | R853C008FG | G853-314-008-3.5-ML |
| 60032265 | G855-313-014-06.0-M, Alpen Diamonds | R855C0145S | G855-313-014-6-MLSS |
| 60032267 | G855-314-010-06.0-M, Alpen Diamonds | R855C010FG | G855-314-010-6-ML |
| 60032269 | G855-314-012-06.0-M, Alpen Diamonds | R855C012FG | G855-314-012-6-ML |
| 60032271 | G855-314-014-06.0-M, Alpen Diamonds | R855C014FG | G855-314-014-6-ML |
| 60032273 | G856-313-016-07.0-M, Alpen Diamonds | R856C0165S | G856-313-016-7-MLSS |
| 60032275 | G856-313-018-08.0-M, Alpen Diamonds | R856C0185S | G856-313-018-8-MLSS |
| 60032277 | G856-314-012-08.0-F, Alpen Diamonds | R856F012FG | G856-314-012-8-F |
| 60032279 | G856-314-012-08.0-M, Alpen Diamonds | R856C012FG | G856-314-012-8-ML |
| 60032281 | G856-314-014-08.0-F, Alpen Diamonds | R856F014FG | G856-314-014-8-F |
| 60032283 | G856-314-014-08.0-M, Alpen Diamonds | R856C014FG | G856-314-014-8-ML |
| 60032285 | G856-314-014-09.0-M, Alpen Diamonds | R856LC014FG | G856-314-014-9-ML |
| 60032287 | G856-314-016-07.0-C, Alpen Diamonds | R856C016FG | G856-314-016-7-ML |
| 60032289 | G856-314-016-07.0-F, Alpen Diamonds | R856F016FG | G856-314-016-7-F |
| 60032293 | G856-314-016-08.0-C, Alpen Diamonds | R856SC016FG | G856-314-016-8-MLX |
| 60032295 | G856-314-016-08.0-F, Alpen Diamonds | R856F016FG | G856-314-016-8-F |
| 60032297 | G856-314-016-08.0-M, Alpen Diamonds | R856C016FG | G856-314-016-8-ML |
| 60032305 | G856-314-018-08.0-C, Alpen Diamonds | R856SC018FG | G856-314-018-8-MLX |
| 60032307 | G856-314-018-08.0-F, Alpen Diamonds | R856F018FG | G856-314-018-8-F |
| 60032309 | G856-314-018-08.0-M, Alpen Diamonds | R856C018FG | G856-314-018-8-ML |
| 60032311 | G856-314-018-09.0-C, Alpen Diamonds | R856LC018FG | G856-314-018-9-MLX |
| 60032313 | G856-314-021-08.0-C, Alpen Diamonds | R856SC021FG | G856-314-021-8-MLX |
| 60032317 | G856-314-021-08.0-M, Alpen Diamonds | R856C021FG | G856-314-021-8-ML |
| 60032323 | G856-314-021-09.0-C, Alpen Diamonds | R856LC021FG | G856-314-021-9-MLX |
| 60032325 | G856-314-021-09.0-F, Alpen Diamonds | R856F021FG | G856-314-021-9-F |
| 60032327 | G856-314-021-09.0-M, Alpen Diamonds | R856C021FG | G856-314-021-9-ML |
| 60032329 | G856-314-025-07.0-C, Alpen Diamonds | R855C025FG | G856-314-025-7-MLX |
| 60032333 | G856-314-025-07.0-M, Alpen Diamonds | R855C025FG | G856-314-025-7-ML |
| 60032335 | G856-314-025-08.0-C, Alpen Diamonds | R856SC025FG | G856-314-025-8-MLX |
| 60032337 | G856-314-025-08.0-M, Alpen Diamonds | R856C025FG | G856-314-025-8-ML |
| 60032341 | G856L-314-014-10.0-M, Alpen Diamonds | R856L10C014FG | G856L-314-014-10-ML |
| 60032343 | G856L-314-016-10.0-C, Alpen Diamonds | R856LC016FG | G856-314-016-9-MLX |
| 60032345 | G856L-314-016-10.0-M, Alpen Diamonds | R856C016FG | G856-314-016-9-ML |
| 60032349 | G856L-314-018-10.0-M, Alpen Diamonds | R856LC018FG | G856-314-018-9-ML |
| 60032357 | G856LT-314-016-10.0-C, Alpen Diamonds | R856LTC016FG | G856L-314-016-10-TSX |
| 60032359 | G856LT-314-018-11.0-C, Alpen Diamonds | R856LTC018FG | G856L-314-018-11-TSX |
| 60032361 | G856N-314-014-08.0-M, Alpen Diamonds | R856NC014FG | G856N-314-014-8-ML |
| 60032367 | G856T-314-014-08.0-M, Alpen Diamonds | R856TC014FG | G856-314-014-8-TSX |
| 60032369 | G856T-314-018-08.0-C, Alpen Diamonds | R856TC018FG | G856-314-018-8-TSX |
| 60032371 | G856T-314-021-09.0-C, Alpen Diamonds | R856TC021FG | G856-314-021-9-TSX |
| 60032373 | G856T-314-025-07.0-C, Alpen Diamonds | R856TC025FG | G856-314-025-7-TSX |
| 60032379 | G856Z-314-018-08.0-M, Alpen Diamonds | ZREX856M018FG | ZREX856-314-018-8-M |
| 60032381 | G858-313-012-08.0-M, Alpen Diamonds | R858C0125S | G858-313-012-8-MLSS |
| 60032383 | G858-314-012-08.0-EF, Alpen Diamonds | R858EF012FG | G858-314-012-8-F |
| 60032385 | G858-314-012-08.0-F, Alpen Diamonds | R858F012FG | G858-314-012-8-F |
| 60032387 | G858-314-012-08.0-M, Alpen Diamonds | R858C012FG | G858-314-012-8-ML |
| 60032389 | G858-314-014-08.0-EF, Alpen Diamonds | R858EF014FG | G858-314-014-8-F |
| 60032391 | G858-314-014-08.0-F, Alpen Diamonds | R858F014FG | G858-314-014-8-F |
| 60032393 | G858-314-014-08.0-M, Alpen Diamonds | R858C014FG | G858-314-014-8-ML |
| 60032395 | G858-314-016-08.0-M, Alpen Diamonds | R858C016FG | G858-314-016-8-ML |
| 60032397 | G859-314-010-10.0-EF, Alpen Diamonds | R859EF010FG | G859-314-010-10-F |
| 60032399 | G859-314-010-10.0-F, Alpen Diamonds | R859F010FG | G859-314-010-10-F |
| 60032403 | G859-314-010-10.0-M, Alpen Diamonds | R859C010FG | G859-314-010-10-ML |
| 60032409 | G859-314-014-10.0-EF, Alpen Diamonds | R859EF014FG | G859-314-014-10-F |
| 60032411 | G859-314-014-10.0-F, Alpen Diamonds | R859F014FG | G859-314-014-10-F |
| 60032415 | G859-314-014-10.0-M, Alpen Diamonds | R859C014FG | G859-314-014-10-ML |
| 60032419 | G859-314-018-10.0-EF, Alpen Diamonds | R859EF018FG | G859-314-018-10-F |
| 60032421 | G859-314-018-10.0-F, Alpen Diamonds | R859F018FG | G859-314-018-10-F |
| 60032423 | G859-314-018-10.0-M, Alpen Diamonds | R859C018FG | G859-314-018-10-ML |
| 60032425 | G860-314-010-04.0-M, Alpen Diamonds | R860C010FG | G860-314-010-4-ML |
| 60032427 | G860-314-012-05.0-EF, Alpen Diamonds | R860SF012FG | G860-314-012-5-F |
| 60032429 | G860-314-012-05.0-F, Alpen Diamonds | R860F012FG | G860-314-012-5-F |
| 60032431 | G860-314-012-05.0-M, Alpen Diamonds | R860SC012FG | G860-314-012-5-ML |

| REF new | Description new | REF old | Description old |
|----------|---------------------------------------|---------------|----------------------|
| 60032433 | G860-314-012-06.5-F, Alpen Diamonds | R860SF012FG | G860-314-012-6.5-F |
| 60032435 | G860-314-012-06.5-M, Alpen Diamonds | R860SC012FG | G860-314-012-6.5-ML |
| 60032439 | G860-314-014-05.0-F, Alpen Diamonds | R860F014FG | G860-314-014-5-F |
| 60032441 | G860-314-014-05.0-M, Alpen Diamonds | R860C014FG | G860-314-014-5-ML |
| 60032443 | G862-313-010-08.0-M, Alpen Diamonds | R862C0105S | G862-313-010-8-MLSS |
| 60032445 | G862-314-010-08.0-EF, Alpen Diamonds | R862EF010FG | G862-314-010-8-F |
| 60032447 | G862-314-010-08.0-F, Alpen Diamonds | R862F010FG | G862-314-010-8-F |
| 60032449 | G862-314-010-08.0-M, Alpen Diamonds | R862C010FG | G862-314-010-8-ML |
| 60032451 | G862-314-012-08.0-EF, Alpen Diamonds | R862EF012FG | G862-314-012-8-F |
| 60032453 | G862-314-012-08.0-F, Alpen Diamonds | R862F012FG | G862-314-012-8-F |
| 60032455 | G862-314-012-08.0-M, Alpen Diamonds | R862C012FG | G862-314-012-8-ML |
| 60032457 | G862-314-014-08.0-F, Alpen Diamonds | R862F014FG | G862-314-014-8-F |
| 60032459 | G862-314-014-08.0-M, Alpen Diamonds | R862C014FG | G862-314-014-8-ML |
| 60032461 | G862-314-016-08.0-EF, Alpen Diamonds | R862EF016FG | G862-314-016-8-F |
| 60032463 | G862-314-016-08.0-F, Alpen Diamonds | R862F016FG | G862-314-016-8-F |
| 60032465 | G862-314-016-08.0-M, Alpen Diamonds | R862C016FG | G862-314-016-8-ML |
| 60032467 | G862Z-314-014-08.0-M, Alpen Diamonds | ZREX862M014FG | ZR862-314-014-8-M |
| 60032469 | G863-314-012-10.0-EF, Alpen Diamonds | R863EF012FG | G863-314-012-10-F |
| 60032471 | G863-314-012-10.0-M, Alpen Diamonds | R863C012FG | G863-314-012-10-ML |
| 60032473 | G863-314-012-10.0-M, Alpen Diamonds | R863C012FG | G863-314-012-10-ML |
| 60032475 | G863-314-016-10.0-EF, Alpen Diamonds | R863EF016FG | G863-314-016-10-F |
| 60032477 | G863-314-016-10.0-F, Alpen Diamonds | R863F016FG | G863-314-016-10-F |
| 60032479 | G863-314-016-10.0-M, Alpen Diamonds | R863C016FG | G863-314-016-10-ML |
| 60032481 | G863Z-314-014-10.0-M, Alpen Diamonds | ZR863ZM014FG | ZR863-314-014-10-M |
| 60032483 | G864-314-014-12.0-F, Alpen Diamonds | R864F014FG | G864-314-014-12-F |
| 60032485 | G864-314-014-12.0-M, Alpen Diamonds | R864C014FG | G864-314-014-12-ML |
| 60032487 | G877-314-010-06.0-M, Alpen Diamonds | R877C010FG | G877-314-010-6-ML |
| 60032489 | G877-314-012-06.0-M, Alpen Diamonds | R877C012FG | G877-314-012-6-ML |
| 60032491 | G877K-314-012-06.0-EF, Alpen Diamonds | R877KC012FG | G877K-314-012-6-ML |
| 60032493 | G877K-314-014-06.0-M, Alpen Diamonds | R877KC014FG | G877K-314-014-6-ML |
| 60032497 | G877K-314-016-06.0-M, Alpen Diamonds | R877KC016FG | G877K-314-016-6-ML |
| 60032499 | G877K-314-018-06.0-M, Alpen Diamonds | R877KC018FG | G877K-314-018-6-ML |
| 60032501 | G878-314-010-08.0-F, Alpen Diamonds | R878F010FG | G878-314-010-8-F |
| 60032503 | G878-314-010-08.0-M, Alpen Diamonds | R878C010FG | G878-314-010-8-ML |
| 60032505 | G878-314-012-08.0-F, Alpen Diamonds | R878F012FG | G878-314-012-8-F |
| 60032507 | G878-314-012-08.0-M, Alpen Diamonds | R878C012FG | G878-314-012-8-ML |
| 60032509 | G878-314-014-08.0-M, Alpen Diamonds | R878C014FG | G878-314-014-8-ML |
| 60032511 | G878-314-016-08.0-M, Alpen Diamonds | R878C016FG | G878-314-016-8-ML |
| 60032513 | G878K-314-010-08.0-F, Alpen Diamonds | R878KC010FG | G878K-314-010-8-F |
| 60032515 | G878K-314-010-08.0-M, Alpen Diamonds | R878KC010FG | G878K-314-010-8-ML |
| 60032517 | G878K-314-012-08.0-F, Alpen Diamonds | R878KC012FG | G878K-314-012-8-F |
| 60032519 | G878K-314-012-08.0-M, Alpen Diamonds | R878KC012FG | G878K-314-012-8-ML |
| 60032521 | G878K-314-014-08.0-F, Alpen Diamonds | R878KC014FG | G878K-314-014-8-F |
| 60032523 | G878K-314-014-08.0-M, Alpen Diamonds | R878KC014FG | G878K-314-014-8-ML |
| 60032525 | G878K-314-016-08.0-C, Alpen Diamonds | R878KC016FG | G878K-314-016-8-MLX |
| 60032529 | G878K-314-016-08.0-M, Alpen Diamonds | R878KC016FG | G878K-314-016-8-ML |
| 60032531 | G878K-314-018-08.0-C, Alpen Diamonds | R878KC018FG | G878K-314-018-8-MLX |
| 60032533 | G878K-314-018-08.0-F, Alpen Diamonds | R878KC018FG | G878K-314-018-8-F |
| 60032535 | G878K-314-021-08.0-C, Alpen Diamonds | R878KC021FG | G878K-314-021-8-MLX |
| 60032537 | G878K-314-021-08.0-M, Alpen Diamonds | R878KC021FG | G878K-314-021-8-ML |
| 60032549 | G878KT-314-016-08.0-C, Alpen Diamonds | R878KTC016FG | G878K-314-016-8-TSX |
| 60032553 | G879-314-012-10.0-M, Alpen Diamonds | R879C012FG | G879-314-012-10-ML |
| 60032555 | G879-314-014-10.0-C, Alpen Diamonds | R879C014FG | G879-314-014-10-ML |
| 60032561 | G879-314-016-10.0-M, Alpen Diamonds | R879C016FG | G879-314-016-10-ML |
| 60032563 | G879K-314-012-10.0-M, Alpen Diamonds | R879KC012FG | G879K-314-012-10-ML |
| 60032565 | G879K-314-014-10.0-M, Alpen Diamonds | R879KC014FG | G879K-314-014-10-ML |
| 60032567 | G879K-314-016-10.0-C, Alpen Diamonds | R879KSC016FG | G879K-314-016-10-MLX |
| 60032569 | G879K-314-016-10.0-M, Alpen Diamonds | R879KSC016FG | G879K-314-016-10-ML |
| 60032571 | G879K-314-018-10.0-C, Alpen Diamonds | R879KSC018FG | G879K-314-018-10-MLX |
| 60032573 | G879K-314-018-10.0-M, Alpen Diamonds | R879KSC018FG | G879K-314-018-10-ML |
| 60032575 | G879K-314-021-10.0-C, Alpen Diamonds | R879KSC021FG | G879K-314-021-10-MLX |
| 60032579 | G879K-314-021-10.0-M, Alpen Diamonds | R879KSC021FG | G879K-314-021-10-ML |
| 60032581 | G879K-314-023-10.0-M, Alpen Diamonds | R879KSC023FG | G879K-314-023-10-ML |
| 60032583 | G879KT-314-014-10.0-M, Alpen Diamonds | R879KTC014FG | G879K-314-014-10-TSX |
| 60032585 | G879KT-314-018-10.0-C, Alpen Diamonds | R879KTC018FG | G879K-314-018-10-TSX |
| 60032589 | G880-314-012-06.0-M, Alpen Diamonds | R880C012FG | G880-314-012-6-ML |
| 60032591 | G880-314-014-06.0-M, Alpen Diamonds | R880C014FG | G880-314-014-6-ML |
| 60032593 | G881-314-010-08.0-M, Alpen Diamonds | R881C010FG | G881-314-010-8-ML |
| 60032601 | G881-314-012-08.0-M, Alpen Diamonds | R881C012FG | G881-314-012-8-ML |
| 60032603 | G881-314-014-08.0-C, Alpen Diamonds | R881SC014FG | G881-314-014-8-MLX |
| 60032607 | G881-314-014-08.0-M, Alpen Diamonds | R881C014FG | G881-314-014-8-ML |
| 60032609 | G881-314-016-08.0-C, Alpen Diamonds | R881SC016FG | G881-314-016-8-MLX |
| 60032611 | G881-314-016-08.0-F, Alpen Diamonds | R881F016FG | G881-314-016-8-F |
| 60032613 | G881-314-016-08.0-M, Alpen Diamonds | R881C016FG | G881-314-016-8-ML |
| 60032615 | G881T-314-014-08.0-M, Alpen Diamonds | R881TC014FG | G881-314-014-8-TSX |
| 60032617 | G881T-314-016-08.0-C, Alpen Diamonds | R881TC016FG | G881-314-016-8-TSX |
| 60032619 | G881Z-314-014-08.0-M, Alpen Diamonds | ZREX881M014FG | ZR881-314-014-8-M |
| 60032621 | G881Z-314-016-08.0-C, Alpen Diamonds | ZREX881C016FG | ZR881-314-016-8-C |

| REF new | Description new | REF old | Description old |
|----------|-----------------|---------|-----------------|
| 60032627 | G882 | | |