

# POLISHING INSTRUMENTS

FOR THE DENTAL PRACTICE





1 COMPREPOL ULTRA				
ISO Ø 1/10 mm	030	050	060	100
Head L in mm	7.5	10	7.5	2.5
Polisher type	2101RA	2103RA	2111RA	2127RA
RA 10 pcs	230055AA	230059AA	230063AA	230069AA
RA 50 pcs		230060AA	230064AA	

2 COMPOSHINE ULTRA				
	030	050	060	100
	7.5	10	7.5	2.5
	2201RA	2203RA	2211RA	2227RA
	230072AA	230076AA	230080AA	230086AA
		230077AA	230081AA	

Composite, Compomere



1 COMPREPOL PLUS						
ISO Ø 1/10 mm	030	040	060	100	110	140
Head L in mm	7.5	10	7.5	2.5	1.6	1.6
Polisher type	2301RA	2303RA	2311RA	2327RA	235G11RA	235G14RA
RA 5 pcs	230291AA	230293AA	230295AA	230297AA	60019877	60019878
RA 10 pcs		230294AA	230296AA			

2 COMPOSHINE PLUS						
	030	040	060	100	110	140
	7.5	10	7.5	2.5	1.6	1.6
	2401RA	2403RA	2411RA	2427RA	245G11RA	245G14RA
	230299AA	230301AA	230303AA	230305AA	60019879	60019880
		230302AA	230304AA			

Composite, Compomere



1 PREBROWN			
ISO Ø 1/10 mm	040	045	060
Head L in mm	7.5	10	9.5
Polisher type	1101RA	1103RA	1111RA
RA 10 pcs	230001AA	230005AA	230013AA
RA 50 pcs		230006AA	

2 HIGREEN			
	040	045	060
	7.5	10	9.5
	1201RA	1203RA	1211RA
	230028AA	230032AA	230040AA

Precious Metal Alloy, Amalgam



1 ZIRCOPOL PLUS			
ISO Ø 1/10 mm	040	060	140
Head L in mm	10	7.5	1.6
Polisher type	9703RA	9711RA	975G14RA
RA 5 pcs	230379AA	230381AA	60022004
RA 10 pcs			

2 ZIRCOSHINE PLUS			
	040	060	140
	10	7.5	1.6
	9803RA	9811RA	985G14RA
	230391AA	230393AA	60022005

CAD/CAM Restoration, Zirconia



1 CERAPRESHINE		
ISO Ø 1/10 mm	050	070
Head L in mm	10	10
Polisher type	3201RA	3211RA
RA 5 pcs		
RA 10 pcs	230114AA	230122AA
Ceramic		

2 CERASHINE		
ISO Ø 1/10 mm	050	070
Head L in mm	10	10
Polisher type	3301RA	3311RA
RA 5 pcs		
RA 10 pcs	230139AA	230147AA
Ceramic		

1 CERAFIN PLUS				
ISO Ø 1/10 mm	040	060	110	140
Head L in mm	10	7.5	2	1.6
Polisher type	9303RA	9311RA	9331RA	935G14RA
RA 5 pcs				
RA 10 pcs	230317AA	230319AA	230325AA	60019881
Ceramic				



2 CERAPRESHINE PLUS					
ISO Ø 1/10 mm	030	040	060	110	140
Head L in mm	7.5	10	7.5	2	1.6
Polisher type	9401RA	9403RA	9411RA	9431RA	945G14RA
RA 5 pcs					
RA 10 pcs	230331AA	230335AA	230337AA	230343AA	60019882
Ceramic					

3 CERASHINE PLUS					
ISO Ø 1/10 mm	030	040	060	110	140
Head L in mm	7.5	10	7.5	2	1.6
Polisher type	9501RA	9503RA	9511RA	9531RA	955G14RA
RA 5 pcs					
RA 10 pcs	230349AA	230353AA	230355AA	230361AA	60019883
Ceramic					



BRUSH			
ISO Ø 1/10 mm	060	040	040
Head L in mm	5.5	6.5	6.5
Polisher type	7001RA	7002RA	7003RA
RA 5 pcs	232002AA	232004AA	232006AA
RA 10 pcs	232001AA	232003AA	232005AA
Ceramic, Composite, Metal			

DIAStrip	
Length x Width	80 mm x 2.5 mm
Polisher type	600UM
5 pcs	232014AA
10 pcs	232015AA
Composite, Compomere, Glass Ionomer Cement, Amalgam	




## Order information

ISO number  
(material, shank, shape, grit size, variable diameter)

Head diameter in 1/10 mm  
Head length  
Polisher description  
Packaging sizes available

ISO 658 204 243 523 ...



Illustration

REF in respective packaging size

ISO Ø 1/10 mm	050
Head L in mm	10
Polisher type	2103RA
RA 10 pcs	230059AA
RA 50 pcs	230060AA

## Head Shapes

ISO Head Shape	Description	ISO Head Shape	Description	Polishing System	RPM	Polishing System	RPM
014	Small Wheel	243 040 10	Point	Comprepol Ultra	5'000–10'000	Cerapreshine	5'000–10'000
030	Cup	243 045-055 010	Flame	Composhine Ultra	5'000–10'000	Cerashine	5'000–10'000
042	Wheel	273	Dome	Comprepol Plus	3'000–8'000	Cerafin Plus	5'000–10'000
044	Wheel	274	Bullet	Composhine Plus	3'000–8'000	Cerapreshine Plus	5'000–10'000
107	Cylinder	292	Torpedo	Prebrown	5'000–10'000	Cerashine Plus	5'000–10'000
114	Cylinder	303	Large Lens	Higreen	5'000–10'000	Brush	5'000–10'000
173	Taper	304	Small Lens	Zircopol Plus	5'000–10'000		
243 030-040 7.5	Minipoint	372	Disc	Zircoshine Plus	5'000–10'000		
		543	ShapeGuard			<b>Specialty Head Shape</b>	<b>RPM</b>
						ShapeGuard	10'000–12'000

## Recommended Speed Ranges

© COLTENE – www.coltene.com



**CH** Coltene/Whaledent AG  
Feldwiesenstrasse 20  
9450 Altstätten / Switzerland  
T +41 71 757 53 00  
F +41 71 757 53 01  
info.ch@coltene.com

**UK** Coltene/Whaledent Ltd.  
The President Suite  
Kendal House, Victoria Way  
Burgess Hill, West Sussex  
RH15 9NF / UK  
T +44 1444 235486  
F +44 1444 870640  
info.uk@coltene.com



**COLTENE**