


# MicroMega One RECI USP



## ➤➤➤ MINIMALLY INVASIVE

1mm wire diameter<sup>(3)</sup>

## ➤➤➤ CUTTING EFFICIENCY & CENTERING ACCURACY

Patented asymmetric cross section and inactive tip combined with C.Wire heat treatment<sup>(2)</sup>

## ➤➤➤ FLEXIBILITY

Respect of the canal anatomy thanks to C.Wire heat treatment<sup>(2)</sup>

## ➤➤➤ SAFETY

Good resistance to cyclic fatigue<sup>(1)</sup>