



**PRESIDENT**<sup>®</sup>  
Impression Tray System


**If they weren't transparent you'd think  
they were made of steel!**

# PRESIDENT®

## Impression Tray System






**Rigid by design**



**Anatomically correct**



**Built-in retention**



Feature	Benefit
Rigid by design	A patented system of intersecting ribs and rails provides unprecedented rigidity in a plastic tray
Anatomically correct	A statistical database accumulated over 20 years has enabled Coltène / Whaledent to define the ideal form for each tray system
Built-in retention	The President Tray System's design aggressively retains all types of impression material.
Transparent	Easy and precise positioning of the tray.
Modifiable form	The tray can individually be adjusted with heat.
Single use	No risk of contamination, no time-consuming cleaning.
Economical	Good value for money

### Product Ordering Information: 1 Bag of 12 Trays

Description	Dimensions		
<b>6901 1 Large Upper</b>	<b>A</b> 83.1 mm/3.270"	<b>B</b> 63.2 mm/2.490"	<b>C</b> 18.7 mm/.738"
<b>6902 2 Large Lower</b>	<b>A</b> 79.9 mm/3.147"	<b>B</b> 62.2 mm/2.450"	<b>C</b> 20.0 mm/.788"
<b>6903 3 Medium Upper</b>	<b>A</b> 75.4 mm/2.968"	<b>B</b> 58.3 mm/2.295"	<b>C</b> 17.3 mm/.683"
<b>6904 4 Medium Lower</b>	<b>A</b> 69.5 mm/2.738"	<b>B</b> 58.5 mm/2.303"	<b>C</b> 19.3 mm/.759"
<b>6905 5 Small Upper</b>	<b>A</b> 69.6 mm/2.714"	<b>B</b> 51.7 mm/2.036"	<b>C</b> 17.1 mm/.671"
<b>6906 6 Small Lower</b>	<b>A</b> 61.7 mm/2.431"	<b>B</b> 50.5 mm/1.990"	<b>C</b> 18.3 mm/.720"
<b>6907 7 Quadrant UL/LR</b>	<b>A</b> 47.1 mm/1.854"	<b>B</b> 58.5 mm/2.303"	<b>C</b> 19.3 mm/.759"
<b>6908 8 Quadrant UR/LL</b>	<b>A</b> 47.1 mm/1.854"	<b>B</b> 58.5 mm/2.303"	<b>C</b> 19.3 mm/.759"
<b>6909 9 Anterior</b>	<b>A</b> 65.3 mm/2.569"	<b>B</b> 39.4 mm/1.554"	<b>C</b> 19.2 mm/.756"

© 2013 Coltène/Whaledent AG – www.coltene.com

#### Coltène/Whaledent AG

Feldwiesenstrasse 20  
9450 Altstätten / Switzerland  
Tel +41 71 757 5300  
Fax +41 71 757 5301  
info.ch@coltene.com

#### Coltène/Whaledent Inc.

235 Ascot Parkway  
Cuyahoga Falls, Ohio 44223 / USA  
Tel +1 330 916 8800  
Fax +1 330 916 7077  
info.us@coltene.com

#### Coltène/Whaledent Ltd.

The President Suite  
Kendal House, Victoria Way  
Burgess Hill, West Sussex  
RH15 9NF / UK  
Tel +44 1444 235486  
Fax +44 1444 870640  
info.uk@coltene.com