



PINDEX® **SPLIT SLEEVE**

**For stable dies that fit just right;
not too loose, not too tight...**

- The unique Pindex sleeve design automatically adjusts to the individual diameter of each pin
- Enhances die stability and consistency of fit

PINDEX®...the world's pinned model standard for more than 40 years.

PINDEX® **SPLIT SLEEVE**

The unique PINDEX® sleeve design automatically adjusts to the individual diameter of each pin. Dramatically reduces the possibility of dies that are too tight or too loose.



Self-adjusts to the diameter of the pin.



No technique modification is necessary.



It provides stable fit with the PINDEX® pin.

PINDEX® Pins & Split Sleeves

PX114S Long Pins & Split Sleeves 1,000 pcs.



PX112S Medium Pins & Split Sleeves 1,000 pcs.



PX 117 Short Pins & Plastic Sleeves 100 pcs.



PX145 Self-Articulating Pins & Plastic Sleeves 100 pcs.



PX145S Plastic Sleeves for Articulating Pins 1,000 pcs.

PINDEX® Accessories

PX119 Carbide Drills 3 pcs.

PX120 Reamer 1 pc.

PX125 Laser Retro-Fit Kit 1 pc.

PX146 Flexible Rubber Base Molds (Full arch, right and left quadrants) 3 pc.

PX147 Flexible Rubber Base Molds (Full arch) 1pc.

PX148 Flexible Rubber Base Molds (Right / left quadrants) 2pcs.

PX200 PINDEX® System Handsaw Complete Kit
 PINDEX® Handsaw 1 pc.
 Saw Blades, .007" (.18mm) 12 pcs.
 Saw Blades, .010" (.25mm) 12 pcs.

PX201 Saw Blades, .007" (.18mm) 12 pcs.

PX301 Saw Blades, .007" (.18mm) 100 pcs.

PX202 Saw Blades, .010" (.25mm) 12 pcs.

PX302 Saw Blades, .010" (.25mm) 100 pcs.

P38 Proxima® Saw Blades 12 pcs.

PINDEX® Laser Drilling Machines

Laser Unit	Voltage
PX113L	115V
P113LCE	220/230V
P113LUK	220/230V



Coltène/Whaledent Inc.
 235 Ascot Parkway
 Cuyahoga Falls, OH 44223 / USA
 Tel. USA & Canada + 1 800 221 3046
 + 1 330 916 8800
 Fax +1 330 645 8704
 info.us@coltene.com

Coltène/Whaledent AG
 Feldwiesenstrasse 20
 9450 Altstätten / Switzerland
 Tel. +41 (0)71 757 53 00
 Fax +41 (0)71 757 53 01
 info.ch@coltene.com

Coltène/Whaledent GmbH + Co. KG **EC REP**
 Raiffeisenstrasse 30
 89129 Langenau / Germany
 Tel. +49 (0)7345 805 0
 Fax +49 (0)7345 805 201
 info.de@coltene.com

www.coltene.com

COLTENE

Made in USA

P/N 30592C