

FILE SEQUENCE – STEP BY STEP CARD

HyFlex™ CM

✓ New file



✓ After use



✓ After autoclaving



✗ Twisted (discard)



Coltène/Whaledent GmbH + Co. KG 

Raiffeisenstraße 30

89129 Langenau/Germany

Coltène/Whaledent Inc.

235 Ascot Parkway

Cuyahoga Falls, OH 44223/USA

www.coltene.com

Use a Slow Speed Handpiece at 500 rpm. Recommended Torque Setting up to 2.5 N-cm (25 mN-m, 250 g-cm)

