

Imagine

whipped cream...



...then realize that it is

JET BLUE BITE

Bite Registration Material

- Stays exactly where you want it
- No resistance to closure, no rebound!
- Rock hard for accuracy, yet flexible for easy trimming
- Pick your ideal speed
 - **Fast** - 30 second working time
40 second oral setting time
 - **Superfast** - 20 second working time
30 second oral setting time

Exceptional Handling with firm, accurate, fast set

JET BLUE BITE fast / superfast bite registration material flows with a smooth whipped-cream consistency exactly where you want it in a predictable quantity and ideal ribbon pattern. It offers a fast oral setting time

and sets rock hard for greater accuracy but is still flexible enough for easy trimming or cutting. The superior handling characteristics are the key for a perfect occlusal fit on the first try.

Direct application of JET BLUE BITE fast



Application of JET BLUE BITE fast



Direct application with the unique Spreader Tip



Ideal working-time



Bite registration



Final bite registration



Easy cutting



Easy trimming



Important and helpful accessories: Unique spreader tip for predictable and consistent dispensing ribbon width and depth

JET BLUE BITE	Working time	Oral setting time	Endhardness	Dimensional stability in %
JET BLUE BITE fast	max. 30 s	min. 40 s	89 Shore A	-0,08%
JET BLUE BITE superfast	max. 20 s	min. 30 s	89 Shore A	-0,12%

Internal measurement

JET BLUE BITE fast 50 ml Single Pack	6491
Cartridge	1 x 50 ml
Mixing Tip Pink	6 pcs.
Spreader Tips	12 pcs.
JET BLUE BITE fast 50 ml Refill	6495
Cartridges	4 x 50 ml
Mixing Tip Pink	12 pcs.
Spreader Tips	12 pcs.

JET BLUE BITE superfast 50 ml Single Pack	6493
Cartridge	1 x 50 ml
Mixing Tip Pink	6 pcs.
Spreader Tips	12 pcs.
JET BLUE BITE superfast 50 ml Refill	6497
Cartridges	4 x 50 ml
Mixing Tip Pink	12 pcs.
Spreader Tips	12 pcs.

JET BLUE BITE Accessories		
50 ml System		
Mixing Tip Pink	40 pcs.	6551
Spreader Tips	100 pcs.	6223
Dispenser 50/75	1 pc.	6900



© 2012 Coltène/Whaledent AG – www.coltene.com

Coltène/Whaledent Ltd.
The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
Tel +44 1444 235486
Fax +44 1444 870640
info.uk@coltene.com

Coltène/Whaledent AG
Feldwiesenstr. 20
9450 Altstätten/Switzerland
Tel. + 41 71 757 53 00
Fax + 41 71 757 53 01
info.ch@coltene.com

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223/USA
Tel. + 1 330 916 8800
Fax + 1 330 916 7077
info.us@coltene.com

 **COLTENE®**