





HyFlex EDM OGSF Sequence

protocol



400 rpm
Max torque
2.5 Ncm



300 rpm
Max torque
1.8 Ncm



400 rpm
Max torque
2.5 Ncm



400 rpm
Max torque
2.5 Ncm

Coltène/Whaledent GmbH + Co. KG
Raiffeisenstraße 30
DE-89129 Langenau
T 00 800 11 0000 12
service@coltene.com

Coltène/Whaledent AG
Feldwiesenstrasse 20
CH-9450 Altstätten
T 00 800 11 0000 12
service@coltene.com

 **COLTENE**



Your Endo Guide

Workshop-Kit

80000383 0723

 **COLTENE**