




Surgitip-endo

Dries Canals

FAST!

roeko®
SurgiTip®-endo




Endodontic aspirator tip used to quickly dry root canals

- Unique, multi-directional tip design offers flexibility without kinking
- Simply place into saliva ejector, insert into canal, and aspirate
- After use, only 1-2 paper points are needed to dry canal
- Supplied sterile, ready to use

 **COLTENE**

roeko[®] SurgiT[®]Tip[®]-endo



	Tip opening \varnothing mm	End opening \varnothing mm	Contents	REF
Surgitip-endo	0.35 mm	6.5 mm	20 pc.	60018280

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, OH 44223 / USA
Tel. USA & Canada + 1 800 221 3046
+ 1 330 916 8800
Fax +1 330 645 8704
info@coltenewhaledent.com

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info@coltenewhaledent.ch

Coltène/Whaledent GmbH + Co. KG **EC REP**
Raiffeissenstrasse 30
89129 Langenau / Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info@coltenewhaledent.de

 **COLTENE**