



Pindex®

Model and Die Fabrication System

Uncompromising model and die precision

Laser Drilling Machine

- Assures parallel channels for quick, easy, and reliable pin placement
- More accurate drill positioning for better die registration
- Visible beam in all lighting for reduced eye fatigue
- Focus-free laser outlasts conventional bulbs
- Manual drilling to pre-set depth

System Accessories

- Complete range of pins and mating sleeves
- Durable, high quality carbide drill
- Flexible rubber base molds
- Handsaw with blades



The benchmark in dental model making

 **COLTENE**

Pindex® Model and Die Fabrication System

Drilling Machine	Technical Data	Order Information	Cat. No.
Overall Dimensions:	13.5 x 7.0 inches (34.3 x 17.8 cm)	USA/Canada/Mexico 115V Europe 230V United Kingdom 230V	PX113L P113LCE P113LUK
Net Weight:	8.25 lbs / 3.7 kg		
Nominal Power:	100 watts		
Motor 115V 60HZ, 3000 rpms Class 1 Laser Product IEC 60825-1-2001		Includes manual drilling unit and accessories: carbide drill, chuck tool, shaft lock tool, reamer, rubber mold set, pin assortment.	

Motor 230V 50HZ, 2500 rpms
Class II Laser Product EN60825-1

For optimum performance of your Pindex System, we recommend the use of genuine Pindex Pins and Accessories. Look for the unique "W" on the end of all brass Pindex Pins.

Pins and Sleeves

PX114S	Long Pins & Split Sleeves	1,000 pcs.	PX145	Self-Articulating Pins & Plastic Sleeves	100 pcs.
PX112S	Medium Pins & Split Sleeves	1,000 pcs.	PX144		1,000 pcs.
PX117	Short Pins & Plastic Sleeves	100 pcs.	PX145S	Plastic Sleeves for Articulating Pins	1,000 pcs.
PX115		1,000 pcs.			

Other Accessories

PX119	Carbide Drills	3 pcs.	PX200	INDEX® Handsaw System Complete Kit	
PX120	Reamer	1 pc.		Pindex® Handsaw	1 pc.
PX124	Bulbs (Mark II)	10 pcs.		Saw Blades, .007" (.18mm)	12 pcs.
PX125	Laser Retro-Fit Kit (For PX113 model only)	1 pc.		Saw Blades, .010" (.25mm)	12 pcs.
PX146	Flexible Rubber Base Molds (Full arch, Right & Left Quadrants)	3 pcs.	PX201	Saw Blades, .007" (.18mm)	12 pcs.
PX147	Flexible Rubber Base Molds (Full arch)	1 pc.	PX301	Saw Blades, .007" (.18mm)	100 pcs.
PX148	Flexible Rubber Base Molds (Right/Left Quadrants)	2 pcs.	PX202	Saw Blades, .010" (.25mm)	12 pcs.
			PX302	Saw Blades, .010" (.25mm)	100 pcs.
			P38	Proxima® Saw Blades	12 pcs.



Coltene/Whaledent Inc. 
235 Ascot Parkway
Cuyahoga Falls, OH 44223 / USA
Tel. USA & Canada + 1 800 221 3046
+ 1 330 916 8800
Fax +1 330 645 8704
info.us@coltene.com
Made in USA

Coltene/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
Tel. +41 (0)71 757 53 00
Fax +41 (0)71 757 53 01
info.ch@coltene.com

Coltene/Whaledent GmbH + Co. KG 
Raiffeisenstrasse 30
89129 Langenau / Germany
Tel. +49 (0)7345 805 0
Fax +49 (0)7345 805 201
info.de@coltene.com

coltene.com

 **COLTENE**

P/N 30591C