


MicroMega HandFiles⁽¹⁾




Pre-sterilized instruments
Blisters of 6 HandFiles

(1) K-Files, K-Reamers, H-Files
micromega.coltene.com

COLTENE's commitment to patient, practitioner and staff safety at the office: pre-sterilized, ready-to-use instruments



- Guaranteed sterilization in compliance with regulatory standards
- Simplified logistics with ready-to-use instruments



The same high product quality, new packaging

- Blister of 6 pre-sterilized instruments
- Multiple-use files



MicroMega K-Files



MicroMega K-Reamers



MicroMega H-Files

K-Files, K-Reamers, H-Files

- Available lengths :
 - > 21 mm
 - > 25 mm
 - > 31 mm **NEW SIZE**
- Size from 006 to 140⁽¹⁾
- Taper 2%



For more information about MicroMega HandFiles, please visit our website : micro-mega.com/catheterism



(1) Scan the QR Code to download the detailed list of references or find it on channel.micro-mega.com

©COLTENE

Legal manufacturer:
Micro-Mega SA
12 rue du Tunnel
25000 BESANCON / France
T +33 3 81 54 42 42
F +33 3 81 54 42 30
commercial@micro-mega.com

MicroMega K-Files - MicroMega K-Reamers - MicroMega H-Files
Advertising | Year of CE-marking: 2020 | Medical device class
Is according to Directive 93/42/EEC | Notified body: LNE/G-
MED | Medical device for dental care, meant for professional
dental use only. See product labelling and instructions for
use on the website: micro-mega.com

80000152-C-01.21

