

RETREATMENT

HyFlex™ Remover

NEW



The essential piece for your retreatment

remover.coltene.com

 **COLTENE**

- One file to remove the root filling material
- Respects the root canal anatomy
- Minimally invasive 1mm wire
- Efficient without any solvent
- Non-cutting tip



ORDER INFORMATION



30/.07

Maximum torque: **2.5 N.cm**

Speed of rotation: **400-800 rpm**

HyFlex Remover, 5 pcs

30/.07 Remover File

Active Blade



19 mm


23 mm

60023649

60023648

15.5 mm

18.5 mm



◀ For best results complete the procedure with HyFlex EDM shaping files. You are also free to use a different file system such as HyFlex CM or system of your choice.



For more information about HyFlex, please take a look at the brochure or visit our website hyflex.coltene.com.

Coltène/Whaledent Inc.
235 Ascot Parkway
Cuyahoga Falls, Ohio 44223 / USA
T +1 330 916 8800
F +1 330 645 8704
info.us@coltene.com

Coltène/Whaledent AG
Feldwiesenstrasse 20
9450 Altstätten / Switzerland
T +41 71 757 5300
F +41 71 757 5301
info.ch@coltene.com

© COLTENE – www.coltene.com

COLTENE

PN: 31624A