

HySolate SyntX Dam

The most advanced latex free dental dam
for faster isolation and confident retraction.




Benefits

- ✓ Extremely elastic and resistant ensuring stability and impeccable isolation
- ✓ Same retraction properties as the best latex dams
- ✓ Resistant to endodontic rinsing solutions: NaOCl (3% and 6%), EDTA and CHX
- ✓ Easy to position even in complex operative situations
- ✓ Blue color enhances contrast for improved visibility
- ✓ Pre-printed on one side
- ✓ 100% latex free, non-allergenic polyisoprene dam

Why HySolate SyntX?

- 1 No tearing during placement
- 2 Reliable isolation every time
- 3 Pre-printed for faster punching



With pre-printed
template

HySolate SyntX Dam is made of **polyisoprene** which combines exceptional elasticity and tear resistance without the allergy risks associated with natural latex.



Clinical case provided by Dr. Alessio Zanza



"I had the pleasure of testing the new HySolate SyntX Dam from COLTENE in advance, and I was truly impressed. Its ease of use, combined with an extremely elastic and resistant material, makes every application quick and precise. Contact points are easily overcome without risk of tearing, while the color optimizes visual focus and is perfect for clinical photography. Moreover, the fact that it is latex-free makes it ideal even for allergic patients. Highly recommended."


Dr. Alessio Zanza, Rome, Italy

Check-out the full HySolate portfolio including clamps, instruments and accessories

Order Information

SIZE	THICKNESS	COLOUR	ARTICLE NO.	SHEETS
152 × 152 mm 6" × 6"	medium	blue	60037567	30
			60037568	75
127 × 127 mm 5" × 5"			60037566	30



We help you preserve and improve natural teeth. The Swiss way. 

 **COLTENE**