

BRILLIANT Crios

Reinforced CAD/CAM composite
for permanent Restoration



Benefits

- ✓ Tooth-like flexibility and controlled wear make it gentler on the antagonist
- ✓ Thermally cured, silicate glass-reinforced composite matrix
- ✓ Minimal post-processing: No firing, no glazing, ultra-quick polishing
- ✓ Easily adjustable and repairable with universal composites of BRILLIANT line

NEW – Xtra Bleach shade for modern aesthetics

Why BRILLIANT Crios?

- 1 Mimics the natural behaviour of tooth structure thanks to shock-absorbing effects
- 2 Protects both the restored tooth and the antagonist
- 3 Highly efficient workflow: Easy chairside milling, minimal post-processing



Shades and Translucencies

Low Translucent



NEW

High Translucent



Super Translucent



„The mechanical properties are great, but what amazes me about BRILLIANT Crios is its aesthetic performance, as it maintains particular brightness in the mouth and can be polished easily.“

Dr. Raphael Monte Alto, Brazil



Order Information

Size	Low (LT)		High (HT)		Super (ST)
	12	14	12	14	14
Xtra Bleach		60042907			
Bleach	60019998	60019985			60024174
A1	60019999	60019986	60020007	60019994	
A2	60020000	60019987	60020008	60019995	
A3	60020001	60019988	60020009	60019996	
A3.5	60020002	60019989			
B1	60020003	60019990	60020010	60019997	
B2	60020004	60019991			
B3	60020005	60019992			
C2	60020006	60019993			
Universal					60024175

More Information



BRILLIANT Crios is also available as Disc.

We help you preserve and improve natural teeth. The Swiss way. 

 **COLTENE**