

KENDA Deluxe

THE LUXURY OF COMPOSITE POLISHING

Developed in close collaboration with clinical experts, KENDA Deluxe is a high-performance, two-step diamond polishing system designed for modern aesthetic dentistry.



More information



Application Video



 COLTENE



One Pre-Polisher – Adaptive & Gentle Smoothing without changing the established surface morphology



One High-Gloss Polisher – Brilliant Results in Seconds



ShapeGuard – One Shape for All Surfaces

ALL COMPOSITES, COMPOMERS & HYBRID - CERAMICS POLISHERS

 DIAMOND POLISHERS



Polisher type 6, 12, 20, 25, 100 pcs	.006		.006		.006		.006		.006		.006		.006							
	.012		.012		.012		.012		.012		.012		.012							
	.025		.025		.020		.025		.025		.020		.025							
	.100		.100		.100		.100		.100		.100		.100							
Ø mm	6.5		3.3		14		5		10		6.5		3.3		14		5		10	
Head L in mm	10		7.4		2		10		5.5		10		7.4		2		10		5.5	
RPM	Rec		7'500 - 10'000		5'000 - 7'000		7'500 - 10'000		7'500 - 10'000		5'000 - 7'500		6'000 - 8'000		5'000 - 7'500		5'000 - 7'500		5'000 - 7'500	
	Max				10'000								10'000							

REF (Set & Samples)

0400.012	1	2	1	1	1	1	2	1	1	1
0400.003	-	-	1	1	-	-	-	1	-	-

Contouring - Rough Surface

Finishing & Polishing - Smooth Surface

Ultra Shine Finish - Glossy Surface



1st Step



2nd Step