

SELF-ADHESIVE DUAL-CURING RESIN CEMENT



SoloCem

One step reliable cementation

**THE NEW
SUPERKRAFT**

 **COLTENE**

Your partner for virtually all indirect bonding indications. Depending on the situation, you can cement self-adhesively or adhesively. For universal use – that's where **SoloCem** comes in.

Advantages at a glance



Only one cement for virtually all indications



Long-lasting restorations



Aesthetic luting



Easy to use



A partner you can count on – always

SoloCem



A host of different indications, yet always the same solution

You believe in speed and efficiency? Then choose SoloCem, because self-adhesive application is not only fast but also safe. However, depending on the clinical case, a stronger adhesive bond may be necessary. Even then, SoloCem is the right choice – use ONE COAT 7 UNIVERSAL as "adhesion booster" in addition.


Less administration.

More time for important things

As you only need one cement for both luting options (adhesive and self-adhesive) with SoloCem, you not only save space in your inventory, but also time spent for the administration of additional products. That means: More time for what is most important – the patient.

Avoids marginal gaps due to low shrinkage

Shrinkage according to Archimedes



Source: internal data


* Not registered trademarks of the COLTENE Group

If shrinkage of a luting cement is high, this may result in marginal gaps. These pose the risk of bacterial infiltration. SoloCem has low shrinkage and is therefore a safe option.

Optimal integration = aesthetic restoration

An aesthetic restoration must meet multiple requirements. On the one hand, the colour of the cement should be inconspicuous and virtually match the natural shade of the teeth. Good colour stability and resistance to foods with a strong discolouration effect - such as red wine or curry - is also important for perfect integration. SoloCem offers both. And owing to the fluorescent properties of SoloCem, fillings are invisible under special lighting conditions, such as black light in a nightclub.


Average discolouration after storage in red wine, curry and cress solution as well as distilled water



Source: Liebermann et al. (2017): The Effect of Different Storage Media on Color Stability of Self-Adhesive Composite Resin Cements for up to One Year

* Not registered trademarks of the COLTENE Group

A simple and straightforward range, but always with the right shade




**THE NEW
SUPERKRAFT**



Strong individually, unbeatable together

Shear bond strength on different materials



High bonding strengths are indicative of a good bond between cement and the various materials. This is crucial for a successful treatment. SoloCem reliably exhibits high bonding values on both dentine and enamel, especially in combination with ONE COAT 7 UNIVERSAL. In self-adhesive application, the material also proves convincing on various other materials such as metals or ceramics.



"The universal self-adhesive resin-based luting cement SoloCem is a help to clinicians."

Prof. Brian Millar, BDS, PHD, FDSRCS, FHEA, FDT, UK

"SoloCem offers the option of achieving a safe and stable long-term result in a simple manner."

Dr. med. dent. Anton Wetzel, Switzerland



Order Information

60013986

SoloCem Intro Kit

- 1 × SoloCem Dentin, 5 ml
- 1 × SoloCem Trans, 5 ml
- 10 × Mixing Tip Brown Short Super Fine
- 10 × Mixing Tip Brown Short Fine

60014062

SoloCem Dentin Refill

- 1 × SoloCem Dentin, 5 ml
- 5 × Mixing Tip Brown Short Super Fine
- 5 × Mixing Tip Brown Short Fine

60014061

SoloCem White Opaque Refill

- 1 × SoloCem White Opaque, 5 ml
- 5 × Mixing Tip Brown Short Super Fine
- 5 × Mixing Tip Brown Short Fine

60019840

SoloCem Trans Refill

- 1 × SoloCem Trans, 5 ml
- 5 × Mixing Tip Brown Short Super Fine
- 5 × Mixing Tip Brown Short Fine

6747

Mixing Tip Brown Short Super Fine

- 40 × Mixing Tip, Ø 1 mm

6759

Mixing Tip Brown Short Fine

- 40 × Mixing Tip, Ø 1.8 mm

© COLTENE – www.coltene.com

Coltène/Whaledent Ltd.

The President Suite
Kendal House, Victoria Way
Burgess Hill, West Sussex
RH15 9NF / UK
T +44 1444 235486
info.uk@coltene.com

Coltène / Whaledent AG

Feldwiesenstrasse 20
9450 Altstätten / Switzerland
T +41 71 757 53 00
info.ch@coltene.com

 COLTENE